Français

..... 4

Deutsch

..... 36

Italiano

..... 68

English

..... 100

Español

..... 132

Nederlands

..... 164

SOMMAIRE

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INTRODUCTION 103

WARM-UP CATEGORY 104

Preparation for seasonal endurance activity (eg: hiking, bicycle touring)	105
Supplement to the Muscle tone, Muscle firming, Shaping or Refinement programmes of the Body sculpt category	105
Achieving peak form before competition for resistance sports (eg: 800 metres, 1500 metres, mountain biking, mountain stage in cycling, cycling pursuit on track, 200 metres swimming, slalom)	105
Prevention of contractures in the muscles of the back of the neck for cyclists	105
Muscular and circulation conditioning before physical activity	106
Use of Potentiation program to optimise the effects of explosive strength (sprints, jumps, squash, football, basketball, etc.) immediately before specific training or competition	106

FITNESS CATEGORY 107

Preparation for a fitness enthusiast who wants to develop harmonious shoulder muscles with a moderate increase in muscle volume	108
Muscle volume for a body-builder	108
Muscle volume for a body-builder training three times a week	109
Muscle volume for a body-builder training five times a week	109
Developing and sculpting your abdominal belt	110
Developing aerobic capacity for an expert fitness enthusiast	110
Preparation for a fitness enthusiast who wants to develop his/her capacity to sustain an intense effort	111
Preparation for a fitness enthusiast who wants to develop and strengthen the large pectoral muscles	111
Preparation for a well trained person who wants to improve the tone of his/her quadriceps using the Step programme	112

CROSS-TRAINING CATEGORY 113

Resuming activity after a break in training: varied muscle preparation	114
Preparation for a fitness enthusiast who wants to adopt a varied muscle training approach	114
Preparation for a long-distance runner who also wants to improve his speed at the end of the race	115
Improving the user's capacity for intense efforts without increasing the risks of injury	115
Competition preparation for an amateur practising a team sport	116

BODY SCULPT CATEGORY 117

To firm and sculpt the body	118
To reduce your waist and tone up your abdominal belt	118
To improve the firmness and shape of your buttocks	118

RECOVERY CATEGORY 119

Practical rules of use for the Undulated massage 1 and Undulated massage 2 programs	121
You want to recover better and faster after an effort	122
Treatment of uncomfortable muscular tensions in the back of the neck	122
Treatment of a localised sensation of heaviness or an occasional state of fatigue	122
Returning to a state of well-being after stressful situations	122
To eliminate unpleasant body sensations and recover an optimal state of relaxation	123
To restore the feeling of muscle equilibrium and soothe undesirable manifestations	123
Use of the Regeneration program to eliminate muscle fatigue more quickly (cross-country running, football, basketball, tennis, etc.) and restore good muscular sensations more quickly	123
Prevention of the feeling of heaviness in the legs	124
Preparation for a fitness enthusiast who wants to optimise the effects of his/her active stretching sessions for quadriceps	124

PAIN RELIEF CATEGORY **125**

Neuralgia of the upper limb (brachial neuralgia)	127
Chronic muscular pain (polymyalgia)	127
Contracture (eg: localized contracture in external side of the calf)	127
Chronic muscular pain in the back of the neck (cervical pain)	127
Muscular pain in the thoracic region (thoracic back pain)	128
Muscular pain in the low back region (low back pain)	128
Sharp and recent muscular pains affecting a muscle in the low back (lumbago)	128
Chronic elbow pain (epicondylitis = tennis elbow)	129
Torticollis (recent and acute back of the neck muscle pains)	129

TEST CATEGORY **130**

Locating a motor point (eg: locating the motor points of the vastus internus and the vastus externus of the quadriceps)	130
--	-----

INTRODUCTION

The applications that follow are given by way of example. They provide a better understanding of how electrostimulation sessions can be combined with voluntary activity. These protocols will help you to determine the best procedure to follow depending on your needs (choice of the program, muscular group, duration, placement of the electrodes, body position).


You can of course choose a muscular group different from the one indicated, depending on the activities you practice or the body area you wish to stimulate.

The specific applications provide you information on the placement of electrodes and the stimulation position to adopt. This information is represented by numbers (placement of electrodes) and letters (stimulation position). They refer to the pictures for the placements of electrodes and the pictograms for the body positions appearing on the flap of the manual's cover.

When the Compex muscle stimulator **is not equipped** with the **mi-SENSOR**, you must select in the device the muscular group indicated on the figurine. When no figurine appears (→), it is not necessary to choose the muscular group: the selection is automatic.

When the Compex muscle stimulator **is equipped** with the **mi-SENSOR**, selection of the muscle group is automatic. Consequently, it is unnecessary to refer to the small pictogram.

Example :

Programs	Effects	Uses	Stimulation energies	mi functions
Capillarization	Very marked increase in blood flow Development of capillaries	In the period before competition, for those participating in endurance or resistance sports To improve endurance in those whose fitness level is not high	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE
Toning massage	Activation of blood circulation Recovery of muscular contractile properties Invigorating effect	To prepare the muscles in an ideal manner before an unusual/one-time physical activity	Progressively increase the stimulation energies until pronounced muscle twitches are obtained Make sure that the stimulation energy is sufficient in order to impose significant muscle contractions	mi-SCAN mi-RANGE
Potentiation	Increased speed of contraction and power gain Less nervous effort to attain maximum strength	For optimum muscle preparation immediately before competition	Maximum tolerable energy (0-999) The energy is adjusted during the phases of the program when the muscle relaxes and reacts with simple muscle twitches	mi-SCAN

Specific applications

Preparation for seasonal endurance activity (eg: hiking, bicycle touring)

Sporting physical activities that last a long time (walking, cycling, cross-country skiing, etc.) require strong muscles with a good capillary circulation so that the muscle fibres can be properly oxygenated. When a long-duration physical activity is not practised, or not practised regularly enough, the muscles lose their capacity to effectively consume oxygen and the capillary system becomes rarefied. This lack of muscular quality makes exercising uncomfortable and limits the speed of recuperation, and is responsible for many unpleasant side effects such as drowsiness, contractures, cramp and swelling.

To restore to muscles their endurance and ensure full activation of the capillary system, the Compex offers a very efficient mode of stimulation.


Cycle duration: 6-8 weeks, before starting the hike, 3 x/week, alternating muscular groups

Program: **Capillarization** 8  and **Capillarization** 25 

Supplement to the Toning, Firming or Shaping programs of the Body sculpt category

Voluntary muscle work or work carried out with electrostimulation can be supplemented to good effect by the **Capillarization** program. Indeed, the marked increase in blood flow in the muscle capillaries produced by the **Capillarization** treatment allows maximum oxygenation of the muscles and surrounding tissue. Muscle fibres are thus correctly oxygenated and obtain maximum benefit from the exercise they have just carried out. In addition, drainage and oxygenation of the surrounding tissue restores equilibrium and rounds off the work.

Cycle duration: 13 weeks, Follow each session **Toning, Firming or Shaping** with 1 **Capillarization** session

Program: **Capillarization**  on the same muscle group as the programs of the Body sculpt category

Achieving peak form before competition for resistance sports (eg: 800 metres, 1500 metres, mountain biking, mountain stage in cycling, cycling pursuit on track, 200 metres swimming, slalom)

Sports that require maximum effort lasting between 30 seconds and 5 minutes are the so-called resistance sports. Rapid fibres must work at a capacity close to their maximum and be capable of maintaining this intensive work throughout the event without weakening, in other words the rapid fibres must be resistant.

The **Capillarization** program, which produces a very significant increase in the blood flow in the muscles, brings about development of the intramuscular capillary system (capillarization). This growth of the capillaries works best if it occurs around the rapid fibres. In this way, the latter's exchange surface with the blood increases, allowing an improved supply of glucose, better diffusion of oxygen and quicker evacuation of lactic acid. Capillarization therefore allows the rapid fibres to be more resistant and to maintain their optimum capacity over a longer period. However, prolonged or too frequent use of this program may induce a modification of rapid fibres into slow fibres, thereby risking a decline in performance for strength and speed sports. It is therefore important to follow closely the recommendations below concerning the specific application in order to benefit from the positive effects of this treatment.

Cycle duration: 1 week before competition, 2 x/day, with a 10 minute break between the 2 sessions

Program: **Capillarization** 8 

Prevention of contractures in the muscles of the back of the neck for cyclists

During voluntary training or competitions, some athletes experience the problem of contractures in support muscles (for example, the back of the neck for cyclists). Use of the **Capillarization** program reduces, and may resolve, this problem. The marked increase in blood circulation and development of the capillary system improve the oxygenation of fibres and their exchanges with the blood.

The muscle therefore becomes less prone to contractures.

Cycle duration: 3 weeks, 1 x/day

Program: **Capillarization 15** 

Muscular and circulation conditioning before physical activity

People who regularly practise sports are well aware of the transition, which is often painful, between a rest activity and an occasional physical effort.

The objective of the usual warming-up techniques is to respond to this need to activate progressively the physiological functions involved in one-time physical activity. This physiological need is also very desirable for moderate but unusual physical efforts, as is frequently the case for most of us (treks, biking, jogging, etc.).

The **Toning massage** program offers an ideal benefit in the form of optimal muscular and circulation conditioning before any type of physical effort. It allows us to avoid the uncomfortable sensations – usually experienced during the first few minutes of unusual physical efforts – and to limit the secondary consequences of insufficient preparation (aches, etc.).

Cycle duration: To be used on the muscles most subject to stress in the physical activity involved (in this example: the quadriceps), in the last thirty minutes preceding physical activity

Program: **Toning massage 8** 

Warm-up

Use of Potentiation program to optimise the effects of explosive strength (sprints, jumps, squash, football, basketball, etc.) immediately before specific training or competition

The **Potentiation** program should be used on the key muscles involved in the discipline practised. In this example, the key muscles of the sprinter (quadriceps) will be stimulated. For other disciplines the choice of muscles to be stimulated may be different.

The **Potentiation** program must not take the place of the voluntary warm-up usually performed before the competition. Activation of the cardiovascular system, short accelerations becoming progressively faster, practice starts and stretching will therefore be carried out by the athlete according to his normal routine.

A **Potentiation** session of short duration (approximately 3 minutes) will be applied on the sprinter's quadriceps immediately before the start of his race (or races, in the event of qualifying events). The specific muscular activation of the **Potentiation** program allows the maximum level of performance to be attained in the very first seconds of the race.

Cycle duration: 1 x

Always carry out the session as soon as possible before the start and at all events not more than 10 minutes before the start. After 10 minutes, the phenomenon of potentiation of muscle fibres quickly wears off

Program: **Potentiation 8G** 

Fitness Category

Programs	Effects	Uses	Stimulation energies	Mi functions
Muscle building	Improved muscle trophism Balanced increase of muscular tone and volume	To improve the musculature in general (strength, volume, tone)	Maximum tolerable energy (0-999)	Mi-SCAN Mi-ACTION
Body building	Increased muscle volume Increased muscle fibre diameter Improved muscle resistance	For body-building enthusiasts and users who want to increase their muscle mass	Maximum tolerable energy (0-999)	Mi-SCAN Mi-ACTION
Muscle definition	Increased muscle tone without a significant volume increase To have firm muscles	For those who want very firm and wiry muscles, but without a great muscle volume increase	Maximum tolerable energy (0-999)	Mi-SCAN Mi-ACTION
Aerobic	Improved aerobic endurance Improved oxygen consumption by muscle fibres	For those who do fitness exercises regularly to increase their endurance and improve the harmony of their figure	Maximum tolerable energy (0-999)	Mi-SCAN Mi-ACTION
Anaerobic	Improved muscle capacity for intense and prolonged effort	For those who want to progress or improve their comfort level for activities that require resistance	Maximum tolerable energy (0-999)	Mi-SCAN Mi-ACTION
Power	Improved contractile capacities of the muscular fibres Improved muscular strength with a slight increase in volume	For fitness enthusiasts who want to improve their muscular strength, an important quality for many physical activities	Maximum tolerable energy (0-999)	Mi-SCAN Mi-ACTION
Step	The effects on muscles are identical to those obtained by using a stepping machine Increased muscle tone	To be used on quadriceps only For muscles that are already trained and prepared by active stepping exercises	Maximum tolerable energy (0-999)	Mi-SCAN Mi-ACTION

■ At the end of a cycle, you may either start a new cycle at the next level up or do some maintenance at the rate of 1 session a week at the last level reached.

Specific applications

Preparation for a fitness enthusiast who wants to develop harmonious shoulder muscles with a moderate increase in muscle volume

Most voluntary physical activities, such as jogging and cycling, require the shoulder muscles to do only limited work. That is why it is particularly important to compensate this under-utilisation by associating Compex sessions with your voluntary training schedule. The **Muscle building** program imposes a lot of specific work to the muscles of the upper part of your body and thus ensures the harmonious development of your shoulders, with firm and well shaped muscles. Unlike voluntary exercises involving heavy weights that can be traumatic for your joints and tendons, stimulation with Compex involves no, or very little stress for your joints and tendons.

Cycle duration: 5 weeks 4 x/week, alternating muscular groups

Progression in the levels:

Weeks 1-5: Step to the next level every week

Eg for 1 week

Mon: 45'-1h of voluntary physical activity (jogging, swimming, cycling, fitness activity, etc.), then 1 x **Muscle building** 17H 

Tue: Rest

Wed: 1 x **Muscle building** 18C 

Thu: 45'-1h of voluntary physical activity (jogging, swimming, cycling, fitness activity, etc.)

Fri: 1 x **Muscle building** 17H 

Sat: Rest

Sun: 1 x **Muscle building** 18C 

Program: **Muscle building** 18C  and 17H 

Muscle volume for a body-builder

Despite repeated efforts during their voluntary training, many body-builders encounter difficulties in developing certain muscle groups. The specific stimulation imposed on muscles by the **Body building** program significantly increases the volume of the stimulated muscles. In addition, for a similar session time, the Compex **Body building** program provides a greater volume gain than voluntary training. The additional training imposed by this stimulation program on muscles not sufficiently receptive to traditional training provides a solution for the harmonious development of all muscle groups without recalcitrant areas.

To obtain optimum progress, you are advised to:

- 1) precede the **Body building** sessions with short voluntary training focused on strength; for example 3 series of 5 repetitions at 90% of maximum force;
- 2) carry out a **Capillarization** (of the Warm-up category) session directly after the **Body building** session

Muscle volume for a body-builder training three times a week

In this example, we assume that the body-builder wants to work on his/her biceps, which he/she considers to be his/her weak point. It is of course possible to stimulate other muscles as well. Also this stimulation routine can be applied to several muscular groups at the same time (e.g. biceps then calves).

Cycle duration: 8 weeks, 3 x/week

Progression in the levels:

Week 1: **Body building** level 1

Weeks 2-3: **Body building** level 2

Weeks 4-5: **Body building** level 3

Weeks 6-8: **Body building** level 4

Eg for 1 week

Mon: Rest

Tue: Voluntary training focused on the muscles of the lower limbs, followed by active work on the biceps: 3 series of 5 repetitions at 90% of Fmax., then 1 x **Body building** 20D  followed by 1 x **Capillarization** 20Q 

Wed: Rest

Thu: Voluntary training focused on the muscles of the trunk, followed by active work on the biceps: 3 series of 5 repetitions at 90% of Fmax., then 1 x **Body building** 20D  followed by 1 x **Capillarization** 20Q 

Fri: Rest

Sat: Voluntary training focused on the muscles of the upper limbs, followed by active work on the biceps: 3 series of 5 repetitions at 90% of Fmax., then 1 x **Body building** 20D  followed by 1 x **Capillarization** 20Q 

Sun: Rest

Programs: **Body building** 20D and **Capillarization** 20Q 

Muscle volume for a body-builder training five times a week

In this example, we assume that the body-builder wants to work on his/her calves that he/she considers to be his/her weak point. It is of course possible to stimulate other muscles as well. Also this stimulation routine can be applied to several muscular groups at the same time (e.g. calves then biceps).

Cycle duration: 12 weeks, 5 x/week

Progression in the levels:

Week 1: **Body building** level 1

Weeks 2-3: **Body building** level 2


Weeks 4-5: **Body building** level 3

Weeks 6-8: **Body building** level 4


Weeks 9-12: **Body building** level 5


Eg for 1 week


Mon: Voluntary training focused on the muscles of the lower limbs, followed by active work on the calf muscles: 3 series of 5 repetitions at 90% of Fmax., then **Body building** 4A  followed by 1 x **Capillarization** 4Q 

Tue: Voluntary training focused on the muscles of the lower limbs, followed by active work on the calf muscles: 3 series of 5 repetitions at 90% of Fmax., then 1 x **Body building** 4A  followed by 1 x **Capillarization** 4Q 

Wed: Rest

Thu: Voluntary training focused on the muscles of the lower limbs, followed by active work on the calf muscles: 3 series of 5 repetitions at 90% of Fmax., then 1 x **Body building** 4A  followed by 1 x **Capillarization** 4Q 

Fri: Voluntary training focused on the muscles of the lower limbs, followed by active work on the calf muscles: 3 series of 5 repetitions at 90% of Fmax., then 1 x **Body building** 4A  followed by 1 x **Capillarization** 4Q 

Sat: Voluntary training focused on the muscles of the upper limbs, followed by active work on the biceps: 3 series of 5 repetitions at 90% of Fmax., then 1 x **Body building** 4A  followed by 1 x **Capillarization** 4Q 

Sun: Rest

Programs: **Body building** 4A and **Capillarization** 4Q 

Developing and sculpting your abdominal belt

To restore or develop the qualities of your abdominal muscles, they must do a lot of intensive work. Also, the voluntary exercises of the abdominal belt required to achieve results are unpleasant and often also dangerous. Indeed, if they are not performed correctly, they can cause crushing of the lumbar spine accompanied, as a consequence, by the onset or aggravation of pain in the lower back.

Complex offers a specific stimulation mode to make your abdominal belt firmer, improve its tone and give it a more harmonious shape, or even to develop a "rock-hard" abdomen. All this can be achieved without the minimum stress and danger for the lumbar spine.

Cycle duration: 10 weeks

Weeks 1-5: 3 x **Muscle building** 10I 
(for the whole abdominal belt) or 11I 
(for the rectus abdominis muscles)/week

Weeks 6-10: 3 x **Muscle definition** 10I 
(for the whole abdominal belt) or 11I 
(for the rectus abdominis muscles)/week

Progression in the levels:

Weeks 1-5: **Muscle building** Step to the next level every week

Weeks 6-10: **Muscle definition** Step to the next level every week

Programs: **Muscle building** and **Muscle definition** 10I or 11I 

Developing aerobic capacity for an expert fitness enthusiast

In this example, the fitness enthusiast wants to work above all on his/her thighs. It is of course possible to stimulate other muscles as well. Also this stimulation routine can be applied to several muscular groups at the same time (e.g. thighs then deltoids).

The intensity of the efforts sustained during many fitness exercises is supported essentially by the aerobic metabolism, for which a good utilisation of oxygen is crucial.

To use oxygen well, it is necessary to have:

- 1) an efficient cardio-vascular system, which can therefore transport large quantities of oxygenated blood to the muscles;
- 2) muscles capable of absorbing the maximum possible amount of the oxygen that reaches them.

The **Aerobic** program enables the user to improve the capacity of the muscles to consume oxygen. It is particularly suited to active and well trained fitness enthusiasts. The latter will therefore improve their comfort level when practising aerobic exercises. Also, the length and/or number of these exercises may thus be gradually increased and progressively done at more sustained rhythms. The effectiveness of the work done will thus greatly increase, with significant muscular benefits.

Cycle duration: 12 weeks, 3 x/week

Progression in the levels:

Week 1: **Aerobic** level 1

Weeks 2-3: **Aerobic** level 2

Weeks 4-5: **Aerobic** level 3

Weeks 6-8: **Aerobic** level 4

Weeks 9-12: **Aerobic** level 5

Eg for 1 week

Mon: 45'-1 h of physical activity in the fitness facility, then 1 x **Aerobic** 8G 

Tue: Rest

Wed: 45'-1 h of physical activity in the fitness facility, then 1 x **Aerobic** 8G 

Thu: Rest

Fri: 1 x **Aerobic** 8G 

Sat: Rest

Sun: 1 h of voluntary physical activity (jogging, swimming, cycling, etc.)

Program: **Aerobic** 8G 

Preparation for a fitness enthusiast who wants to develop his/her capacity to sustain an intense effort

In this example, the fitness enthusiast wants to work above all on his/her quadriceps. It is of course possible to stimulate other muscles as well. Also this stimulation routine can be applied to several muscular groups at the same time (e.g. thighs then buttocks).

Resistance efforts are characterised by high levels of work intensity sustained for as long as possible, often up to the point of exhaustion. During such efforts, it is mainly the lactic anaerobic metabolism that comes under stress. The factor that limits performance is muscular fatigue linked, among other things, to the accumulation of lactic acid and the capacity of the muscle to tolerate a high concentration of this acid.

The specific routine of muscular fibre activation of the **Anaerobic** program produces a major improvement in the tolerance of the muscle for high concentrations of lactic acid. The user thus improves his/her resistance: an effort can be sustained for a longer time at the same level of intensity or for the same time at a higher level of intensity.

Cycle duration: 3-6 weeks, 2-3 x/week

Eg for 1 week

Mon: 1 x **Anaerobic** 8G 

Tue: Normal training in a fitness facility

Wed: 1 x **Anaerobic** 8G 

Thu: Normal training in a fitness facility

Fri: 1 x **Anaerobic** 8G  (optional)

Sat: Normal training in a fitness facility or other outdoor physical activity

Sun: Rest

Program: **Anaerobic** 8G 

Preparation for a fitness enthusiast who wants to develop and strengthen the large pectoral muscles

To obtain well developed pectoral muscles and increase their strength, these muscles must be subjected to intensive work. In voluntary training, this type of work involves exercises with heavy loads in a fitness centre. But this kind of work can be traumatic for the joints and tendons involved.

The development and strengthening of the pectoral muscles can be useful for sports activities requiring efficient upper limbs and shoulders. The **Power** program is the perfect solution for working on these muscles to improve their contractile qualities, thus contributing to the development of their contraction strength. Using this programme, you can also achieve a harmonious development of muscular volume.

Cycle duration: 3-6 weeks, 3 x/week

Program: **Power** 43B 

Preparation for a well trained person who wants to improve the tone of his/her quadriceps using the Step programme

The **Step** program imposes on the quadriceps do the same kind of work as voluntary exercises on a stepping machine. The contractions occur at a very fast rate, which accelerates from level to level, thus progressively increasing the difficulty of the sessions. This kind of very tonic work improves the firmness of the muscle tissue and helps the user to develop well shaped and very resistant thighs.

The use of the **Step** program should be limited only to muscles that are well prepared.

The **Step** program can be used in two ways:

Mode n° 1: The classic approach with progress from one level to the next with cycles lasting several weeks each;

Mode n° 2: During the same session, the five levels are used for a few minutes each in increasing order.

Cycle duration:

Mode n° 1: 3-6 weeks., 3 x/week

Mode n° 2: 6 weeks, 3 x/week

Week 1: 1 session = 6 minutes with each session of 5 levels

Week 2: 1 session = 7 minutes with each session of 5 levels

Week 3: 1 session = 8 minutes with each session of 5 levels

Week 4: 1 session = 9 minutes with each session of 5 levels

Week 5: 1 session = 10 minutes with each session of 5 levels

Week 6: 1 session = 12 minutes with each session of 5 levels

Program: **Step** 8G –

Cross-training Category

Programs	Effects	Uses	Stimulation energies	Mi functions
Combined workout	Training and preparation of muscles for all kinds of muscular work through different working sequences Maintaining good general physical condition	At the beginning of the season to "re-set" the muscles after a rest period and before more intensive and specific training During the season for the persons who do not have any preferences for a single kind of performance and prefer to exercise their muscles with different workout methods	Maximum tolerable energy (0-999)	Mi-SCAN Mi-ACTION
Cross-endurance	Training and preparation of muscles for all kinds of muscular work through different working sequences, concentrating on work of moderate intensity and long duration	For persons who want, above all, to improve their endurance, but also to develop other muscular qualities Varied training concentrating on the development of aerobic muscular qualities	Maximum tolerable energy (0-999)	Mi-SCAN Mi-ACTION
Cross-resistance	Training and preparation of muscles for all kinds of muscular work through different working sequences, concentrating on high-intensity work practised for as long as possible Increased capacity of the fibres to support work leading to an accumulation of lactic acid	For persons who want to intensify their training schedule without increasing the risks of injury To improve the resistance of muscular fibres to fatigue	Maximum tolerable energy (0-999)	Mi-SCAN Mi-ACTION
Cross-explosivity	Training and preparation of muscles for all kinds of muscular work through different working sequences, with an emphasis on the development of contraction speed Improvement of muscular explosive capacity	For persons who want to subject their muscles to different types of work, while concentrating on the improvement of contraction speed	Maximum tolerable energy (0-999)	Mi-SCAN Mi-ACTION

At the end of a cycle, you may either start a new cycle at the next level up or do some maintenance at the rate of 1 session a week at the last level reached.

Specific applications

Resuming activity after a break in training: varied muscle preparation

In this example, the fitness enthusiast wants to work above all on the latissimus dorsi. It is of course possible to stimulate other muscles as well. Also this stimulation routine can be applied to several muscular groups at the same time (e.g. latissimus dorsi then thighs).

An increasing number of people emphasise a diversified training carried on in a natural setting. During training, different workout routines are alternated and different muscles are stimulated. For example, after a few minutes of slow jogging, the trainee does a series of accelerations followed by some jumps, before resuming a slower jog, and so on. The objective is to work on different muscular qualities without however emphasising any one in particular. This kind of activity is frequently done either at the beginning of the season, for general muscular reactivation, or done regularly by leisure-time sports enthusiasts who do not want to emphasise any particular kind of muscular performance, but rather want to stay fit and reach the right level in all kinds of muscular work.

At the beginning of the season or after a significant break in training, the resumption of physical and/or sports activity should be done progressively and become increasingly specific. It is thus normal to perform a few initial sessions with the aim of making the muscles do all kinds of work to prepare them for later training that will be more intensive and more oriented towards a specific kind of performance.

Through these 8 sequences that automatically follow each other, the **Combined workout** program imposes different kinds of work to the stimulated muscles and thus makes them used to all kinds of effort.

Cycle duration: 3-6 weeks, 3 x/week

Program: **Combined workout** 18C 

Preparation for a fitness enthusiast who wants to adopt a varied muscle training approach

In this example, the fitness enthusiast wants to work above all on her arms. It is of course possible to stimulate other muscles as well. Also this stimulation routine can be applied to several muscular groups at the same time (e.g. arms then thighs).

An increasing number of people emphasise a diversified training carried on in a natural setting. During training, different workout routines are alternated and different muscles are stimulated. For example, after a few minutes of slow jogging, the trainee does a series of accelerations followed by some jumps, before resuming a slower jog, and so on. The objective is to work on different muscular qualities without however emphasising any one in particular. This kind of activity is frequently done either at the beginning of the season, for general muscular reactivation, or done regularly by leisure-time sports enthusiasts who do not want to emphasise any particular kind of muscular performance, but rather want to stay fit and reach the right level in all kinds of muscular work.

Most fitness enthusiasts want to improve their muscle qualities and, with this aim in mind, pursue objectives that may be different: improved muscle tone, increased endurance, greater muscle volume, etc. The means used to achieve these different objectives involve different kinds of voluntary training and the use of specific electrostimulation programs.

Other fitness enthusiasts prefer to develop all their muscular qualities, without any emphasis on any one of them in particular. In this case, it is useful to subject the muscles, during the same session, to various kinds of stimulation, such as those produced by the **Combined workout** program through its 8 routines that automatically form a sequence.

Cycle duration: 3-6 weeks, 2-3 x/week

Eg for 1 week

Mon: 1 x **Combined workout** 21D 

Tue: Normal training in a fitness facility

Wed: 1 x **Combined workout** 21D 

Tu: Normal training in a fitness facility

Fri: 1 x **Combined workout** 21D 

Sat: Normal training in a fitness facility or other outdoor physical activity

Sun: Rest

Program: **Combined workout** 21D 

Preparation for a long-distance runner who also wants to improve his speed at the end of the race

Any kind of physical activity of moderate intensity practised for a long period involves the aerobic metabolism. A good aerobic system depends on an efficient cardiovascular system. At the same time, it requires muscles capable of consuming the oxygen that reaches them (oxidative capacity). Electrostimulation develops this capacity and thus improves the endurance of the stimulated muscular groups.

However, it is useful, during a stimulation session, to change the focus of the program by alternating endurance with other aspects, such as speed, strength and resistance. This allows you not only to develop the endurance of your muscular fibres, but also to work simultaneously on other aspects that may be necessary when practising your chosen physical activity. A long-distance runner, for example, may need strength to cope with the slopes along the route or to increase his speed at the end of the race.

The characteristics of its sequences make the **Cross-endurance** program particularly suitable for improving endurance, but also for improving contraction speed, strength and resistance.

Cycle duration: 12 weeks, 3 x/week

Progression in the levels:

Week 1: **Cross-endurance** level 1

Weeks 2-3: **Cross-endurance** level 2

Weeks 4-6: **Cross-endurance** level 3

Weeks 7-8: **Cross-endurance** level 4

Weeks 9-12: **Cross-endurance** level 5

Eg for 1 week

Mon: Voluntary training: warm up 20', then 1-2 series of 6 x {30" fast / 30" slow}

Slow jogging 10' at the end of the session, then 1 session **Capillarization** 8D 

Tue: 1 session **Cross-endurance** 8G 

Wed: Loose jogging 60', then 1 session **Capillarization** 8D 

Tu: 1 session **Cross-endurance** 8G 

Sa : Loose jogging 45-60'

Di : Extended run 1h30 (moderate speed), then 1 session **Capillarization** 8D 

Programs: **Cross-endurance** 8G and **Capillarization** 8D 

Improving the user's capacity for intense efforts without increasing the risks of injury

In this example, we assume that the person wants to work above all on his calves. It is of course possible to stimulate other muscles as well. Also, the **Cross-resistance** program can be used on several muscle groups at the same time, that is to say the same stimulation procedure can be applied to both calves and buttocks, for instance.

Fitness enthusiasts who regularly do several sessions of intensive voluntary activity every week frequently run the risk of over-stressing the muscles and joints involved. The **Cross-resistance** program offers a good alternative for intensifying the training schedule without increasing the risks of injuries. Also, the alternation of varied stimulation sequences means that contrasting types of work requiring several different muscular qualities can be imposed on the muscles, while also improving their capacity to support better the lactic acid accumulated during high-intensity physical effort.

Cycle duration: 3-6 sem., 3 x/sem.

Eg for 1 week

Mon: 1h of physical activity in the fitness centre

Tue: 1h of physical activity in the fitness centre, then 1 session **Cross-resistance** 4A 

Wed: Rest

Thu: 1h of physical activity in the fitness centre, then 1 session **Cross-resistance** 4A 

Fri: 1 session **Cross-resistance** 4A 

Sat: 1h-1h30 of voluntary physical activity (jogging, swimming, cycling)

Sun: Rest

Program: **Cross-resistance** 4A 

Competition preparation for an amateur practising a team sport

In this example, we assume that the person wants to work above all on his quadriceps. It is of course possible to stimulate other muscles as well. Also, the **Cross-explosivity** program can be used on several muscle groups at the same time, that is to say the same stimulation procedure can be applied to both quadriceps and buttocks, for instance.

Team sports, such as, for example, football, rugby, handball or basket-ball, require a good level of development of different muscular qualities. During a match, strength and explosive strength are decisive qualities, but a good cardio-vascular system plays an important part. Preparation for this type of physical activity should therefore not be concentrated exclusively on a single quality. Rather, it should take into account, through different exercises, the need for diversified voluntary preparation, which also involves factors such as muscular explosive strength, endurance and resistance.

The **Cross-explosivity** program activation routines allow varied work on the muscular groups involved and the selection of different working routines concentrating on muscular explosive strength, an important factor in most team sports.

Cycle duration: 10 weeks, 4 x/week

Progression in the levels:

Weeks 1-2: **Cross-explosivity** level 1

Weeks 3-4: **Cross-explosivity** level 2

Weeks 5-6: **Cross-explosivity** level 3

Weeks 7-8: **Cross-explosivity** level 4

Weeks 9-10: **Cross-explosivity** level 5

Eg for 1 week

Mon: Collective training, then 1 session **Cross-explosivity** 8G 

Cross-training

Tue: 1 session **Cross-explosivity** 8G 

Wed: Training 45'-1h in a fitness centre

Thu: Rest

Fri: Collective training, then 1 session

Cross-explosivity 8G 

Sat: Rest

Sun: 1 session **Cross-explosivity** 8G 

Program: **Cross-explosivity** 8G 

Body sculpt Category

Programs	Effects	Uses	Stimulation energies	mi functions
Toning	To tone the muscles	To be used to begin with in order to tone up and prepare the muscles before more intensive firming work	Maximum tolerable energy (0-999)	mi-SCAN mi-ACTION
Firming	To regain muscle firmness and restore the support function of the muscles	To be used as the main treatment for muscle firming	Maximum tolerable energy (0-999)	mi-SCAN mi-ACTION
Shaping	To define and sculpt the body when the muscles are already firm	To be used when the firming phase is finished	Maximum tolerable energy (0-999)	mi-SCAN mi-ACTION
Abs	Combine the beneficial effects of the Toning and Firming programs adapted to the specificities of the muscles of the abdominal belt	To reduce your waist line, while improving the muscular support of the abdominal wall	Maximum tolerable energy (0-999)	mi-SCAN
Buttocks	Subjects the buttocks muscles to the most appropriate work for improving the figure	To restore and improve the tone of the buttocks region	Maximum tolerable energy (0-999)	mi-SCAN

Specific applications

To firm and sculpt the body

Toning, Firming and Shaping programs impose perfectly adapted and progressive work on the muscles. This very intensive muscular activity (hundreds of seconds of sustained contractions) first tones and firms your muscles, then redefines their contours to sculpt your body.

In this example, we assume that the fitness enthusiast wants to work on her arms. It is of course possible to stimulate other muscles as well. Also this stimulation routine can be applied to several muscular groups at the same time (e.g. arms then thighs).

Cycle duration: 13 weeks

Weeks 1-3: 4 x **Toning** 21D /week

Weeks 4-8: 4 x **Firming** 21D /week

Weeks 9-13: 4 x **Shaping** 21D /week

Progression in the levels:

Weeks 1-3: **Toning** Step to the next level every session

Weeks 4-8: **Firming** Step to the next level every week

Weeks 9-13: **Shaping** Step to the next level every week

Programs: **Toning, Firming and Shaping** 21D 

To reduce your waist and tone up your abdominal belt

The **Abs** program is especially appropriate for improving the qualities of support and firmness of the abdominal belt muscles. The various work routines, associated with the tonic massage performed during the rest phases, guarantee beneficial results on the waist line.

Cycle duration: 3-6 weeks, 4 x/week

Program: **Abs** 49I –

To improve the firmness and shape of your buttocks

The considerable amount of progressive and specific work of the buttocks muscles produced by the **Buttocks** program effectively increases the firmness of these muscles. The shape of the buttocks improves progressively, with beneficial effects on the whole figure.

Cycle duration: 3-6 weeks, 4 x/week

Program: **Buttocks** 9A B or L –

Recovery Category

Programs	Effects	Uses	Stimulation energies	Mi functions
Active recovery	Marked increase in blood flow Accelerates the elimination of waste from the muscular contraction Endorphinic effect (see Pain category) Relaxing effect	To improve and accelerate the muscle recovery after an intensive exercise To be used during the first 3 hours after each session of intensive training or after a competition	Increase the stimulation energies progressively until marked muscle twitches are obtained	Mi-SCAN Mi-RANGE
Relaxing massage	Decreased muscular tension Drainage of the toxins responsible for the exaggerated increase of muscular tone Effect of well-being and relaxation	To eliminate uncomfortable or painful sensations, following an exaggerated increase of muscular tone	Increase the stimulation energies progressively until marked muscle twitches are obtained	Mi-SCAN Mi-RANGE
Reviving massage	Considerable increase of blood flow in the stimulated region Improvement of tissue oxygenation Elimination of free radicals	To contrast effectively sensations of fatigue and localised heaviness	Increase the stimulation energies progressively until marked muscle twitches are obtained	Mi-SCAN Mi-RANGE
Anti-stress massage	Improvement of tissue vascularisation Decreased muscular tension associated with stress Return to calm	To restore a sensation of well-being after stressful situations	Increase the stimulation energies progressively until marked muscle twitches are obtained	Mi-SCAN Mi-RANGE
Undulated massage 1	Activation of the local circulation and improvement of exchanges at tissue level Reduction of unpleasant sensations	To eliminate uncomfortable body tensions and recover or achieve an optimal state of relaxation	Increase the stimulation energies progressively until marked muscle twitches are obtained	Mi-SCAN Mi-RANGE
Undulated massage 2	Improvement of venous return flow Relief of muscular tensions Sedative effect General relaxation	To recover sensations of well-being and relaxation	Increase the stimulation energy progressively until the most pleasant level is reached, while taking care, if possible, to obtain clearly visible muscular twitches or until a pronounced tickling sensation can be felt under the electrodes	Mi-SCAN Mi-RANGE

Programs	Effects	Uses	Stimulation energies	mi functions
Regeneration	Analgesic effect through the release of endorphins Marked increase in blood flow encouraging oxygenation and drainage Activation of the oxydative metabolism Reactivation of proprioceptive pathways	To be used the day after competition as recovery training or as a complement to this type of training, the intensity of which can then be reduced	Increase the stimulation energies progressively until marked muscle twitches are obtained After 10 minutes, when tetanic contractions begin, the energies should be increased to the bearable maximum (0-999) After this phase of tetanic contractions, reduce the energy but ensure that muscle twitches continue to remain pronounced	mi-SCAN mi-ACTION
Heavy legs	Increased venous return Increased supply of oxygen to muscles Drop in muscular tension Elimination of tendency to have cramp	To remove the feeling of heaviness in the legs arising in unusual situations (standing for long periods, heat, hormonal imbalance linked to the menstrual cycle, etc.)	Increase the stimulation energies progressively until marked muscle twitches are obtained	mi-SCAN mi-RANGE
Stretching	To optimise the effects of the voluntary technique of stretching, by means of a reduction of muscular tonus obtained by specific activation of the antagonist of the stretched muscle (reciprocal inhibition reflex)	For all athletes who wish to maintain or improve their muscular elasticity To be used after all training sessions or during a specific stretching session	Increase progressively the intensities to obtain a tension that is clearly perceived	mi-SCAN

Specific applications

Practical rules of use for the Undulated message 1 and Undulated message 2 programs

The Undulated message 1 and Undulated message 2 programs use a new stimulation mode that produces a particularly pleasant wave-like effect, thus increasing the effectiveness of these programs. In most of the sequences of these two programs, stimulation is alternated on the four channels to reproduce faithfully different manual massage techniques (which must always be done in the direction of the venous return flow, that is to say from the extremities toward the heart). This is especially important to obtain the most benefit from the circulatory effects of these techniques and programs. To take full advantage of this wave-like effect, particular care should be paid to the positioning of the electrodes, and the rules for their use indicated below should be followed.

Rule n° 1: number of stimulation channels to use

The correct use of the Undulated message 1 and Undulated message 2 programs requires the use of the 4 stimulation channels and, therefore, of 8 electrodes.

These 4 stimulation channels must be positioned:

- on a limb or a segment of a limb (e.g., a thigh),
- on the back,
- on two limbs (e.g., the two calves).

Rule n° 2: required order of the stimulation channels

The wave effect must necessarily follow the direction of the venous return flow. That is why the required order of the four stimulation channels must be strictly complied with.

The order of the stimulation channels is as follows:

Channel 1: cable fitted with a transparent (mi-SENSOR) or blue colour indicator.


Channel 2: cable fitted with a green colour indicator,

Channel 3: cable fitted with a yellow colour indicator,


Channel 4: cable fitted with a red colour indicator.

It is always recommended to follow the general indications below:


Application on a limb or a segment of a limb

Application on the back

Application on the two limbs

You want to recover better and faster after an effort

A session with the **Active recovery** program, which accelerates the drainage of toxins, will help you to relax your muscles and to restore your muscular balance faster, whatever type of activity you practice.

The program is usually used between 15 minutes and 3 hours after a physical effort. If the result is insufficient, the session can be repeated immediately a second time. It is also possible to repeat a session the following day, if it is found that aches and/or spasms persist.

Cycle duration: To be used throughout the season

Program: **Active recovery** 8◇ 

Treatment of uncomfortable muscular tensions in the back of the neck

Remaining in a seated position for long periods, associated with repetitive movements of the upper limbs (as is frequently the case in front of a computer screen), is often responsible for an uncomfortable – or even painful – increase of muscle tension of the back of the neck muscles. Any other condition of stress may also create situations leading to a state of excessive muscular tension, which is often responsible for painful or uncomfortable sensations.

The in-depth effects produced by the **Relaxing message** program enable to contrast effectively such painful sensations, with a particularly significant result in terms of relaxation.

Cycle duration: To be used on the back of the neck muscles, whenever an occasional sensation of painful muscular tension is felt; to be repeated, if necessary, in cases of particularly strong muscular tension

Program: **Relaxing message** 15◇ 

Treatment of a localised sensation of heaviness or an occasional state of fatigue

The many stresses of everyday life are often responsible for uncomfortable, or even painful, physical sensations. A circulation slowdown is often the result of insufficient physical activity, frequently aggravated by the need in one's working life to remain in the same position for many hours (in a seated position, for example).

Although not serious, this simple "vascular slowdown" is, however, frequently the cause of unpleasant sensations (for example: a sensation of heaviness, often localised in the lower limbs, but also, at times, in another region of the body).

The **Reviving message** program produces, in the greatest comfort, a reactivation of the blood circulation that allows an acceleration of tissue oxygenation and the elimination of painful sensations, due to insufficient physical effort.

Cycle duration: To be used on the calf muscles, whenever an occasional sensation of heaviness is felt; to be repeated, if necessary, in cases of persistent discomfort

Program: **Reviving message** 25◇ 

Returning to a state of well-being after stressful situations

The many stressful situations of everyday life and the worries they entail create a certain degree of psychic tension to which each of us reacts differently. This hyperactivity of the brain and the emotions also has repercussions on our body, which may respond in many unwanted ways.

Undoubtedly, the most frequent problem is a sensation of discomfort associated with persistent muscular tension. This is sometimes intensified by physical constraints, such as the need to remain for a long time in positions unsuited to the well-being of our body. For a person under stress, it is difficult to relax. The muscles are always too tense, producing uncomfortable and even painful sensations. In this kind of situation, the **Anti-stress message** program can reduce muscular tensions through low-frequency stimulation producing a sensation of perfect calm. Applied to the affected region, at the rate of two sessions a day, this program elimi-

nates the unpleasant sensations produced by a stressful situation and produces a relaxing effect leading to a return to a state of well-being.

Cycle duration: Use this program daily during the stressful period or whenever a painful muscular sensation is felt

Program: **Anti-stress message** 15◇ 

To eliminate unpleasant body sensations and recover an optimal state of relaxation

To relax and eliminate uncomfortable physical sensations are objectives pursued by many people, especially after a difficult and often tiresome day.

To achieve these objectives, the **Undulated message 1** program offers an ideal solution for all those people who want to improve their well-being in this way. The automatic succession of the different program sequences allows the user to obtain beneficial selective effects to overcome the various habitual components of physical discomfort: circulatory slow-down, exaggerated muscular tone, micro-stresses on the joints, etc. Also, the wave-like effect obtained through the alternation of the stimulation phases on the 4 channels further increases the effectiveness of the program and produces a pleasant sensation that is particularly appreciated.

The program can be used effectively on all parts of the body; however, the back with the muscles that support the spinal vertebrae are particularly prone to these problems. The following indication presents a specific application to relieve this region.

Cycle duration: Use this treatment at will, whenever you want to and/or need it

Program: **Undulated message 1** 45◇ 

To restore the feeling of muscle equilibrium and soothe undesirable manifestations

Staying for a long time in uncomfortable static positions and/or continually performing tasks that involve a high level of physical stress, both associated with a state of psychological stress, are frequently responsible for physical "disorders".

Although generally speaking benign, all too frequently they entail psycho-bodily distress in our modern society.

Too much brain-work, usually referred to as "stress", and not enough or too much physical stress lead to this state of fatigue, tension or heaviness that produces a sensation of general discomfort.

The **Undulated message 2** program offers a highly beneficial relaxation session for all those people who are faced occasionally or regularly with these negative situations. The judicious combination of the wave-like effect with the various stimulation sequences, each taking advantage of the recognized benefits of the latest advances in the medical field, allow the achievement of matchless results.

Cycle duration: Use this treatment at will, whenever you want to and/or need it

Program: **Undulated message 2** 46◇ 

Use of the Regeneration program to eliminate muscle fatigue more quickly (cross-country running, football, basketball, tennis, etc.) and restore good muscular sensations more quickly

The **Regeneration** program, known also as the "day-after program" must be used on key muscles for the discipline being practised. In this example, key muscles for the cross-country runner (quadriceps) will be stimulated. For other disciplines, the choice of muscles to be stimulated may be different. Note that this type of program is particularly beneficial for all sports where competitions are repeated at frequent intervals, during tournaments and cups in various sports.

The stimulation session using the **Regeneration** program must be carried out the day after a competition to replace or supplement so-called "restoration" training, which can therefore be less intensive. Contrary to the **Active recovery** program,

which provokes no tetanic contractions and which must be used during the three hours after the competition or intensive training, the **Regeneration** program is a form of light training which, in addition to an analgesic effect and an increased blood flow, aims to impose a small degree of anaerobic training and likewise to provoke slight tetanic contractions that are not tiring, making it possible to reactivate the proprioceptive pathways. Energy pathways are also gently stimulated, allowing their metabolic equilibrium to re-establish itself.

The session is composed of 6 stimulation sequences that follow on automatically:

1st sequence : analgesic effect

2nd sequence : marked increase in blood flow

3rd sequence : tetanic contractions to restore muscular sensation

4th sequence : activation of the oxydative metabolism

5th sequence : marked increase in blood flow

6th sequence : relaxing effect

Cycle duration: To be used throughout the season, according to the frequency of competitions; the day after each competition.

Program: **Regeneration** 8G 

Prevention of the feeling of heaviness in the legs

The feeling of heaviness in the legs is due to a temporary deficiency in the venous return, and not to major organic lesions. Stagnation of the blood mass in the legs is encouraged in certain situations: prolonged periods in a standing position, long periods continuously seated, intense heat, etc. The insufficient supply of oxygen to the tissue, and in particular the muscles, that results from this situation, causes the feeling of heaviness and discomfort in the legs. Use of the **Heavy legs** program accelerates the return of blood to the veins and has an important relaxing effect on the painful muscles.

Cycle duration: Use this treatment whenever an occasional feeling of heaviness in the legs occurs

Program: **Heavy legs** 25 J-

Recovery

Preparation for a fitness enthusiast who wants to optimise the effects of his/her active stretching sessions for quadriceps

The **Stretching** program consists in stimulating the antagonist (i.e. opposite) muscle of the muscle subjected to stretching in order to exploit a well-known physiological mechanism: reciprocal inhibition reflex. This reflex, based on the proprioceptive sensitivity of muscles, consists of a very marked muscular relaxation. This allows more efficient stretching, since it is carried out on a more relaxed muscle.

Stimulation is therefore performed on the muscle opposing the muscle being stretched. This stimulation consists of a progressive appearance and disappearance (slowly and lasting a long time) of contractions, with complete rest between contractions. It is during the contraction (which increases with each new level) that the athlete stretches the chosen muscle group by means of a traditional voluntary stretching technique.

In this example, stimulation is applied to the hamstrings muscles to facilitate the stretching of the quadriceps during the contraction phases.

Cycle duration: Throughout the season, according to the frequency of voluntary stretching sessions

According to the normal duration of stretching (this depends on schools of physical training and each person's feelings), choose the level that seems most appropriate

The proposed duration of stretching is:

10 seconds for level 1

12 seconds for level 2

14 seconds for level 3

16 seconds for level 4

18 seconds for level 5

Program: **Stretching** 5 

The starting position for the active stretching exercise

Pain relief Category

Table of the pathologies

Pathologies	Programs	References
Neuralgia of the upper limb (brachial neuralgia)	Modulated TENS	See application page 127
Chronic muscular pain (polymyalgia)	Endorphinic	See application page 127
Contracture (eg: localized contracture in external side of the calf)	Decontracture	See application page 127
Chronic muscular pain in the back of the neck (cervical pain)	Cervical pain	See application page 127
Muscular pain in the thoracic region (thoracic back pain)	Thoracic back pain	See application page 128
Muscular pain in the low back region (low back pain)	Low back pain	See application page 128
Sharp and recent muscular pains affecting a muscle in the low back (lumbago)	Lumbago	See application page 128
Chronic elbow pain (epicondylitis = tennis elbow)	Epicondylitis	See application page 129
Recent and acute back of the neck muscle pains (torticollis)	Torticollis	See application page 129

Programs	Effects	Uses	Stimulation energies	mi functions
Modulated TENS	Blocks transmission of pain by the nervous system	All acute or chronic localized pain	Progressively increase the stimulation energies until a pronounced tickling sensation can be felt under the electrodes	mi-TENS
Endorphinic	Analgesic action through the release of endorphins Increased blood flow	To combat chronic muscular pain	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE
Decontracture	Reduced muscular tension Relaxing effect	To combat recent and localized muscular pain	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE
Cervical pain	Analgesic action through the release of endorphins Increased arterial flow	Analgesic current specifically adapted to pains in the back of the neck	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE
Thoracic back pain	Analgesic action through the release of endorphins Increased blood flow	Analgesic current specifically adapted to pain in the thoracic region (below the back of the neck and above the low back)	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE
Low back pain	Analgesic action through the release of endorphins Increased blood flow	Analgesic current specifically adapted to persistent low back pain (lumbar region)	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE
Lumbago	Reduced muscular tension Relaxing effect	Analgesic current specifically adapted to sharp and sudden low back pain (lumbar region)	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE
Epicondylitis	Blocks transmission of pain by the nervous system	Analgesic current specifically adapted to persistent pain in the elbow	Progressively increase the stimulation energies until a pronounced tickling sensation can be felt under the electrodes	mi-TENS
Torticollis	Reduced muscular tension Relaxing effect	Analgesic current specifically adapted to acute and sharp pains in the back of the neck	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE

■ The use of **Pain relief** programs should not be prolonged without medical advice. If the pain is great and/or persistent, you are advised to consult a doctor. Only a doctor can make an accurate diagnosis and put in place therapeutic actions designed to encourage the disappearance of the disorder.

Specific applications

Neuralgia of the upper limb (brachial neuralgia)

Some people suffer from arthritis in the joints of the vertebrae at the back of the neck, or from arthritis or periarthritis of the shoulder. These situations often give rise to pain that runs down one arm and is known as “brachial neuralgia”. These arm pains that start in the shoulder or the back of the neck can be reduced with the **Compex Modulated TENS** program by following the practical recommendations outlined below.

Cycle duration: 1 week, 1 x/day minimum, then adapt according to how the pain develops

According to requirements, the **Modulated TENS** program can be repeated a number of times during the same day

Program: **Modulated TENS 35** ⚡

Chronic muscular pain (polymyalgia)

Some people suffer from muscular pain that affects several muscles or parts of different muscles at the same time. The localization of this chronic pain can vary over time. These continual and diffuse muscle pains are the result of chronic contractures in which acids and toxins accumulate, irritating the nerves and causing pain. The **Endorphinic** program is particularly effective against such pain since, in addition to its pain-reducing effect, it increases the blood flow to the contracted muscular groups and removes from them accumulations of acids and toxins.

Example: localized pain in the biceps. However, the pain may affect other muscular groups. The practical application illustrated below remains valid, but it is then necessary to place the electrodes on the muscular group concerned.

Cycle duration: 4 weeks, 2 x/day, with a 10 minute break between the 2 sessions

You are advised to consult your doctor if no improvement is observed after the first week of use

Program: **Endorphinic 20** ⚡

Contracture (eg: localized contracture in external side of the calf)

After tiring muscle work, intense training or a sporting competition, certain muscles or certain parts of muscles often remain tense and slightly painful. These are known as muscular contractures which should disappear after a few days with rest, good rehydration, a balanced food intake with mineral salts and application of the **Decontracture** program. The phenomenon of contracture frequently affects the calf muscles, but can also occur with other muscles. In this case, simply follow the same practical recommendations as below, but place the electrodes on the muscular group concerned.

Cycle duration: 1 week, 1 x/day

You are advised to consult your doctor if no improvement is observed after the first week of use

Program: **Decontracture 24** ⚡

Chronic muscular pain in the back of the neck (cervical pain)

Positions in which the muscles of the back of the neck remain tense for long periods of time, for example working in front of a computer screen, may be responsible for the onset of pain in the back of the neck or on either side of the base of the neck, in the upper back. These pains are due to a contracture of the muscles, of which the prolonged tension crushes the blood vessels and prevents the blood from supplying and oxygenating the muscle fibres. If this phenomenon is prolonged, there is an accumulation of acid and the blood vessels atrophy. The pain then becomes continuous or appears after only a few minutes spent working in an unfavourable position.

These chronic pains in the back of the neck can be treated effectively with the **Cervical pain** program, which reactivates the circulation, drains accumulations of acid, oxygenates the muscles, develops the capillaries and relaxes the contracted muscles.

Cycle duration: 4 weeks, 2 x/day, with a 10 minute break between the 2 sessions

You are advised to consult your doctor if no improvement is observed after the first week of use

Program: **Cervical pain** 15↵

Muscular pain in the thoracic region (thoracic back pain)

Vertebral arthritis and positions in which the muscles of the spinal column remain under tension for long periods of time are often responsible for the onset of pain in the middle of the back, which is accentuated with fatigue. Pressure of the fingers on the muscles on either side of the spinal column can often trigger a sharp pain.

These pains are due to a contracture of the muscles, of which the prolonged tension crushes the blood vessels and prevents the blood from supplying and oxygenating the muscle fibres. If the phenomenon is prolonged, there is an accumulation of acid and the blood vessels atrophy. The pain then becomes continuous or can appear after only a few minutes spent working in an unfavourable position.

These chronic back pains can be effectively treated with the **Thoracic back pain** program which reactivates the circulation, drains accumulations of acid, oxygenates the muscles, develops the capillaries and relaxes the contracted muscles.

Cycle duration: 4 weeks, 2 x/week with a 10 minute break between the 2 sessions

You are advised to consult your doctor if no improvement is observed after the first week of use

Program: **Thoracic back pain** 13↵

Pain relief

Muscular pain in the low back region (low back pain)

Low back pain is the most frequently encountered pain. In a standing position, the entire weight of the trunk is concentrated on the joints between the last vertebrae and the sacrum. The low back region is therefore under particular strain. The discs between the vertebrae are crushed and the low back muscles contractured and painful.

There are very many treatments for the relief of low back pain sufferers; among these, the specific currents of the **Compex Low back pain** program provide an appreciable improvement and can even resolve the problem if it is essentially muscular in origin.

Cycle duration: 4 weeks, 2 x/day, with a 10 minute break between the 2 sessions

You are advised to consult your doctor if no improvement is observed after the first week of use

Program: **Low back pain** 12↵

Sharp and recent muscular pains affecting a muscle in the low back (lumbago)

During back movement, for example when lifting something, when turning or when standing up straight after bending down, a sudden pain can be triggered in the low back. Those who suffer from this problem present a contracture of the low back muscles and feel a sharp pain in this region; because they cannot stand completely straight, they remain bent over on one side. All of these symptoms indicate what is known as lumbago, which is principally the result of a sharp and intense contracture of the low back muscles (lumbar region).

In a situation such as this, it is always necessary to consult a doctor to receive appropriate treatment.

In addition to such treatment, the specific **Compex Lumbago** program can help effectively to relax the muscles and remove the pain.

Cycle duration: 1 week, 1 x/day

You are advised to consult your doctor if no improvement is observed after the first week of use

Program: **Lumbago** 33↵

Pain relief

Chronic elbow pain (epicondylitis = tennis elbow)

All tendons of the muscles that allow us to stretch our hand, wrist and fingers end at the small external bone mass of the elbow (epicondyle). Hand and finger movements therefore transmit tensions that are concentrated in the tendon endings at this bone mass.

When hand movements are repetitive, as is the case for painters, tennis players or even those who constantly use the mouse of a computer system, small injuries, accompanied by inflammation and pain, develop in the region of the epicondyle. This is known as "epicondylitis", which is characterised by pain around the external bone mass of the elbow when pressure is applied or when the forearm muscles are contracted.

The **Compex Epicondylitis** program supplies specific current to combat this type of pain. It acts effectively as a supplement to rest. However, it is necessary to consult your doctor if the pain gets worse or does not disappear quickly after a few sessions.

Cycle duration: 1 week, 2 x/week minimum, then adapt according to how the pain develops

According to requirements, the Epicondylitis program can be repeated several times during the same day

Program: **Epicondylitis** 36↵

Torticollis (recent and acute back of the neck muscle pains)

As a consequence of various circumstances, such as, for example, a rather sudden movement of the neck, an incorrect sleeping position or exposure to a cold environment, a sharp pain may be felt in the back of the neck often accompanied by a considerable reduction of the mobility of the cervical region. A person affected by this problem has contracted muscles in the back of the neck and feels a sharp pain there. Together these symptoms represent what is known as "torticollis", which is mainly the result of a sharp and intense contracture of the muscles in the back of the neck.

In this kind of situation, it is always necessary to consult a doctor in order to receive appropriate treatment.

In addition, **Compex specific Torticollis** program can effectively help to relax the muscles and eliminate the pain.

Cycle duration: 1 week, 1 x/day

You are advised again to consult your doctor if no improvement is observed after the first week of use

Program: **Torticollis** 15↵

Test Category

Programs	Effects	Uses	Stimulation energies	mi functions
Motor point	Stimulation at low frequencies allowing the extremely precise identification of the ideal position for placing the positive electrode (motor point) The precise positioning of the positive electrode on the motor point of the muscle to be stimulated is an important factor for the effectiveness of treatment	Recommended before any initial muscular stimulation session The motor points thus precisely identified may be easily located again if they are marked using a skin marker pencil	Progressively increase the energy, while manually moving the pen tip on the conductive gel, without ever interrupting skin contact A maximum energy between 5 and 15 is usually sufficient to identify the motor point precisely	None
Demo	Refer to the quick start guide, "Test Your Compex in 5 Minutes"			mi-SCAN mi-TENS mi-ACTION mi-RANGE

Specific applications

Locating a motor point (eg: locating the motor points of the vastus internus and the vastus externus of the quadriceps)

Electrostimulation subjects the muscles to work. The progress achieved depends on the kind of work to which the muscles are subjected, that is to say the program chosen. The electrical pulses generated by these programs are transmitted to the muscles (via the motor nerve) through self-adhesive electrodes. The positioning of the electrodes is one of the determinant factors for ensuring a comfortable electrostimulation session.

It is therefore essential to dedicate special care to this aspect. The correct placement of the electrodes and the use of significant energy allow a large number of muscle fibres to work. The greater the energy, the greater the spatial involvement, that is to say the number of fibres working, and therefore the greater the number of fibres that make progress.

A stimulation lead wire consists of two electrodes:

- a positive electrode (+): red connection,
- a negative electrode (-): black connection.

The positive electrode must be fixed precisely on the motor point of the muscle. The motor points correspond to an extremely localised area where the motor nerve is more excitable. Although the location of the various motor points is now well known, there may nevertheless be variations, which can extend to several centimetres, between different individuals.


The **Motor point** program, associated with the use of the motor point pen, allows to determine with great accuracy the exact location of the motor points for each individual and thus ensure the greatest effectiveness of the programs.

It is advisable to use this program before any initial muscular electrostimulation session. The motor points once located can be easily identified using a skin-marker pencil or in any other way, thus avoiding the need to repeat this search before each session.

Advice on use

To locate the motor point of the vastus internus, proceed as described below:

- 1) Apply a large electrode at the top of the thigh (see picture n° 7),
- 2) Connect the negative terminal (black connection) of channel 1 to the output of this large electrode located towards the internal surface of the thigh,
- 3) Spread a thin but uniform layer of Compex conductive gel in the position indicated for the positive electrode positioned over the vastus internus (see picture n° 7), spreading the gel a few extra centimetres in all directions,
- 4) Connect the positive connection (red) to the tip of the motor point pen and bring the tip of the pen in contact with the conductive gel,

- 5) Select the **Motor point** program and set the area for the thighs; then, start the program by pressing the "+" / "-" keys of channel 1,


6) Very progressively increase the energy of channel 1, until a value between 5 and 15 is reached, while continuously moving the pen tip over the gel layer, but without ever losing contact with the skin gel, to avoid triggering off an electrode fault message,

7) As soon as you observe a muscular response in form of twitching, you have located the vastus internus motor point. Visually locate this motor point and apply a small electrode that should be centred over the motor point. If this causes an electrode fault, ignore the message and do not exit the program: continue with the rest of the operations normally,

8) Remove the pen from the positive connection (red) of channel 1 and connect the small electrode to this connection.

To locate the motor point of the vastus externus, proceed as described below:

- 1) Connect the negative terminal (black connection) of channel 2 to the other output of the large electrode that should be left in place,
- 2) Spread a thin but uniform layer of Compex conductive gel in the position indicated for the positive electrode positioned over the vastus externus (see picture n° 7), spreading the gel a few extra centimetres in all directions,
- 3) Connect the positive connection (red) to the tip of the motor point pen and bring the tip of the pen in contact with the conductive,

- 4) Select the **Motor point** program and set the area for the thighs; then, start the program by pressing the "+" / "-" keys of channel 2,

5) Very progressively increase the energy of channel 2, until a value between 5 and 15 is reached, while continuously moving the pen tip over the gel layer, but without ever losing contact with the skin gel, to avoid triggering off an electrode fault message,

6) As soon as you observe a muscular response in form of twitching, you have located the vastus externus motor point. Visually locate this motor point and apply a small electrode that should be centred over the motor point. If this causes an electrode fault, you can exit the program.

7) Remove the pen from the positive connection (red) of channel 2 and connect the small electrode to this connection.

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