# Physio 5

















**CefarCompex** 

	PRACTICAL GUIDE .....	1
	GUIDE PRATIQUE .....	34
	PRAKTISCHER LEITFADEN ...	68
	GUIDA PRATICA .....	102
	GUÍA PRÁCTICA .....	136

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# CONTENTS

<b>I Background .....</b>	<b>3</b>
<b>II General information .....</b>	<b>4</b>
1 - Electrodes .....	4
2 - Stimulation energy - Intensity .....	4
3 - Body position .....	4
<b>III Common Treatment .....</b>	<b>5</b>
1 - Rehabilitation .....	5
2 - Pain relief .....	15
<b>IV Specific Treatment .....</b>	<b>21</b>
1 - Direct current .....	21
A) Precautions prior to Iontophoresis treatment .....	21
B) Preparing the patient and the area to be treated by iontophoresis .....	21
C) Preparing the electrodes and solution of ionized medicine .....	21
D) Attaching the electrodes .....	21
E) During treatment .....	22
F) After treatment .....	22
2 - Denervated .....	23
A) Introduction .....	23
B) Factors guiding the therapeutic approach .....	23
C) Practical therapeutic approach .....	24
3 - Aesthetic .....	25
4 - Incontinence .....	27
<b>V Conditioning .....</b>	<b>31</b>

**Dear Customer!**

In your hand you have CefarCompex's Practical Guide with specific information for every programme in your stimulator.

The Practical Guide covers all the programmes of CefarCompex's Physio 5. This information is easy to read and its purpose is to serve as an extra support, if needed. You'll find the recommended indication and what effect the treatment will have on your patient. Read the guidelines for the recommended intensity level, and see suggested electrodes placements on the fold out cover.

You could always find more information on CefarCompex's website: [www.cefarcompex.com](http://www.cefarcompex.com)

# I Background

## NMES

NMES (NeuroMuscular Electrical Stimulation) is used successfully both in medical rehabilitation as well as a complement to athletic training at all levels.

The goal of electrical muscle stimulation is to achieve contractions or twitches of the muscles. Normal muscular activity is controlled by the central and peripheral nervous systems, which transmit electrical signals to the muscles. NMES works similarly but uses an external source (the stimulator) with electrodes attached to the skin for transmitting electrical impulses into the body. The impulses stimulate the nerves to send signals to a specifically targeted muscle, which reacts by contracting, just as it does with normal muscular activity.

Electrical muscle stimulation is suitable for all the muscles in the body. It can be used to strengthen muscles weakened by surgery, a fracture, etc., and improve mobility. It is also an excellent tool for stroke rehabilitation, helping patients in handgrip and gait training.

## TENS

TENS (Transcutaneous Electrical Nerve Stimulation) produces good results in acute and chronic pain conditions of many kinds. It is clinically proven and used daily by physiotherapists, other caregivers and top athletes around the world.

TENS activates the pain-inhibiting mechanisms of the nervous system. Electrical impulses from electrodes, placed on the skin over or near the painful area, stimulate the nerves to block the pain signals that are sent to the brain, and the pain is less perceived.

TENS is a safe treatment method and has, in contrast to oral painkillers and other pain relief methods, no known side effects. It may be sufficient as the only treatment form, but it is also a valuable complement to other pharmacological and/or physical treatments.

## II General information

### 1) Electrodes

To receive maximum results from your stimulation, it is important to place the electrodes correctly. For muscle contractions, a precise placement on the motor point ensures maximum comfort and efficiency. On the fold outs of this guide you will find comprehensive illustrations of various electrodes placements.

Before placing the electrodes on your body, be sure to always wash and dry the skin in order to ensure good conductivity. The durability of the electrodes depends on how good the care and maintenance instructions are followed. Ensure that your patients always have their own personal electrodes.

Choose electrode size according to the size of the muscle, use small electrodes on smaller muscles and large electrodes on larger muscles. If the stimulation feels unpleasant, this could be a sign that the contact between the skin and the electrode is insufficient. The electrodes could be worn out and need to be replaced.

Note: Always be sure that the stimulation is turned off (energy level 0) before removing electrodes from the skin.

### 2) Stimulation energy - Intensity

For pain relief treatment with a stimulation frequency between 50-100 Hz, the ground rule is that the stimulation should be perceived stronger than the pain itself. Increase the intensity until you obtain a strong tingling sensation that is not painful.

For pain relief with a stimulation frequency between 2-10 Hz, the intensity shall be high enough to give visible muscle response. Note that the patient will feel the stimulation before there is an actual visible contraction!

For NMES and muscle stimulation, the objective for the therapist is to motivate the patient to tolerate his/her highest possible non-painful stimulation energy level. The higher the stimulation energy, the greater percentage of motor units recruited. In some cases, it could be beneficial to combine stimulation with voluntary contractions/co-contractions. Try to get a progression in the energies reached within a session as well as from one session to the next.

### 3) Body position

You will find an illustration of the recommended body position under the electrode placement illustration. The optimum position depends on the type of programme and which body part you wish to stimulate. When you use massage programmes such as Relaxing massage it is important that the patient is in a position where he can relax properly.

## III Common treatment

### 1) Rehabilitation

#### DISUSE ATROPHY

<b>When?</b>	After a long period of immobilization or diminished movement, a normally innervated muscle suffers from a decrease in volume.
<b>Why?</b>	Slow twitch fibres (type I) are particularly affected by this disuse.
<b>How?</b>	By using frequencies creating tetanic contractions of the Type I fibres. In order to modify the programme as the patient progresses, use level 1 for the first two weeks before changing into level 2 the subsequent weeks.
<b>Electrodes</b>	Electrodes positioned as indicated according to the stimulated muscle. See suggested electrodes placements on the fold out cover.
<b>Intensity</b>	Maximum tolerable stimulation energy, which is one of the key factors for the effectiveness of the treatment. The higher the stimulation energy, the greater the number of muscular fibres (motor units) working. Try to gradually increase the energy level during a treatment session.
<b>Option 2+2</b>	If the patient experiences discomfort and pain during stimulation, the 2+2 function allows this programme (ch 1+2) to be combined with a pain relief programme (ch 3+4).

#### DISUSE ATROPHY, LEVEL 1 (25 MIN)

	Warm-up	Contraction	Active rest	Final recovery
Frequency	6 Hz	35 Hz	4 Hz	3 Hz
Duration of ramp-up	1.5 s	1.5 s	0.5 s	1.5 s
Duration of phase	2 min	6 s	7 s	3 min
Duration of ramp-down	2 s	0.75 s	0.5 s	3 s

#### DISUSE ATROPHY, LEVEL 2 (25 MIN)

	Warm-up	Contraction	Active rest	Final recovery
Frequency	6 Hz	45 Hz	4 Hz	3 Hz
Duration of ramp-up	1.5 s	1.5 s	0.5 s	1.5 s
Duration of phase	2 min	6 s	5 s	3 min
Duration of ramp-down	2 s	0.75 s	0.5 s	3 s

## PREVENTION OF DISUSE ATROPHY

<b>When?</b>	After surgery or a bone fracture, a limb or a section thereof is immobilized, the muscles of this part of the body suffer rapidly from disuse atrophy.
<b>Why?</b>	The rapid decrease of muscle volume is mainly due to a reflex inhibition phenomenon or a total absence of any type of exercise. It's important to note that this type of atrophy concentrates on the slow twitch muscle fibres (Type I).
<b>How?</b>	In order to compensate, NMES could reproduce a series of contractions comparable to the normal level of activity during a day. The treatment time of this programme is relatively long and the programme uses frequencies targeting the slow twitch muscle fibres.
<b>Electrodes</b>	Electrodes positioned as indicated according to the stimulated muscle. See suggested electrodes placements on the fold out cover.
<b>Intensity</b>	Maximum tolerable stimulation energy, which is one of the key factors for the effectiveness of the treatment. The higher the stimulation energy, the greater the number of muscular fibres (motor units) working. Try to gradually increase the energy level during a treatment session.
<b>Option 2+2</b>	If the patient experiences discomfort and pain during stimulation, the 2+2 function allows this programme (ch 1+2) to be combined with a pain relief programme (ch 3+4).

### PREVENTION OF DISUSE ATROPHY, LEVEL 1 (54 MIN)

	Warm-up	Contraction	Active rest	Final recovery
Frequency	6 Hz	30 Hz	4 Hz	3 Hz
Duration of ramp-up	1.5 s	3 s	1.5 s	1.5 s
Duration of phase	2 min	5 s	14 s	3 min
Duration of ramp-down	2 s	1.5 s	1.5 s	3 s

### PREVENTION OF DISUSE ATROPHY, LEVEL 2 (47 MIN)

	Warm-up	Contraction	Active rest	Final recovery
Frequency	6 Hz	40 Hz	4 Hz	3 Hz
Duration of ramp-up	1.5 s	3 s	0.5 s	1.5 s
Duration of phase	2 min	6 s	12 s	3 min
Duration of ramp-down	2 s	0.75 s	0.5 s	3 s

## REINFORCEMENT

<b>When?</b>	Muscular reinforcement is used to strengthen healthy muscles (with normal volume). Use on previously disuse atrophied muscles that have regained their volume.
<b>Why?</b>	To obtain increased strength and stability.
<b>How?</b>	By using frequencies creating tetanic contractions of the fast twitch muscle fibres (type II). Often used when rehabilitation is in progress and normal muscular volume is achieved by previous use of disuse atrophy programmes. We recommend the use of level 1 for the first two weeks and level 2 for subsequent weeks.
<b>Example</b>	Strengthen lateral peroneal muscles to prevent sprains. Strengthen shoulder muscles to prevent shoulder dislocations.
<b>Electrodes</b>	Electrodes positioned as indicated according to the stimulated muscle. See suggested electrodes placements on the fold out cover.
<b>Intensity</b>	Maximum tolerable stimulation energy, which is one of the key factors for the effectiveness of the treatment. The higher the stimulation energy, the greater the number of muscular fibres (motor units) working. Try to gradually increase the energy level during a treatment session.
<b>Option 2+2</b>	If the patient experiences discomfort and pain during stimulation, the 2+2 function allows this programme (ch 1+2) to be combined with a pain relief programme (ch 3+4).

REINFORCEMENT , LEVEL 1 (20 MIN)				
	Warm-up	Contraction	Active rest	Final recovery
Frequency	6 Hz	75 Hz	4 Hz	3 Hz
Duration of ramp-up	1.5 s	1.5 s	0.5 s	1.5 s
Duration of phase	2 min	4 s	10 s	3 min
Duration of ramp-down	2 s	0.75 s	0.5 s	3 s

REINFORCEMENT , LEVEL 2 (20 MIN)				
	Warm-up	Contraction	Active rest	Final recovery
Frequency	6 Hz	85 Hz	4 Hz	3 Hz
Duration of ramp-up	1.5 s	1.5 s	0.5 s	1.5 s
Duration of phase	2 min	4 s	8 s	3 min
Duration of ramp-down	2 s	0.75 s	0.5 s	3 s

## NEURO REHAB (SLOW START)

<b>When?</b>	NMES is an excellent complement to traditional physiotherapy for several neurological diagnostics like stroke. The treatment could be used passively, but preferably used actively in the training situation.
<b>Why?</b>	To aid facilitation and relearning of motor skills
<b>How?</b>	This programme has a slow 4 seconds contraction ramp-up time and a long rest phase. Experience shows that this design of the programme is more suitable for patients within this diagnosis group.
<b>Example</b>	Increase function of a subluxated shoulder. Facilitate and increase function in lower arm - hand.
<b>Electrodes</b>	Electrodes positioned as indicated according to the stimulated muscle. See suggested electrodes placements on the fold out cover.
<b>Intensity</b>	Maximum tolerable stimulation energy, which is one of the key factors for the effectiveness of the treatment. The higher the stimulation energy, the greater the number of muscular fibres (motor units) working. Try to gradually increase the energy level during a treatment session.
<b>Option 2+2</b>	If the patient experiences discomfort and pain during stimulation, the 2+2 function allows this programme (ch 1+2) to be combined with a pain relief programme (ch 3+4).

### NEURO REHAB (SLOW START), LEVEL 1 (20 MIN)

	Warm-up	Contraction	Active rest	Final recovery
Frequency	6 Hz	35 Hz	-	3 Hz
Duration of ramp-up	1.5 s	4 s	-	1.5 s
Duration of phase	2 min	5 s	15 s	3 min
Duration of ramp-down	2 s	2 s	-	3 s

### NEURO REHAB (SLOW START), LEVEL 2 (20 MIN)

	Warm-up	Contraction	Active rest	Final recovery
Frequency	6 Hz	45 Hz	-	3 Hz
Duration of ramp-up	1.5 s	4 s	-	1.5 s
Duration of phase	2 min	5 s	15 s	3 min
Duration of ramp-down	2 s	2 s	-	3 s

## BACK-TRUNK/STABILIZATION

<b>When?</b>	For patients with insufficiency in the back and trunk due to long term pain or a neurological disorder. This programme can be combined with active movements.
<b>Why?</b>	To stabilize the trunk and back muscles. Increase the awareness and postural control.
<b>How?</b>	By stimulating the abdominal or lumbar muscle groups.
<b>Example</b>	Abdominal and/or lumbar stimulation to increase the awareness of pelvic motion. Increase stability and function in the lumbar/abdominal muscles.
<b>Electrodes</b>	Electrodes positioned as indicated according to the stimulated muscle. See suggested electrodes placements on the fold out cover.
<b>Intensity</b>	Maximum tolerable stimulation energy, which is one of the key factors for the effectiveness of the treatment. The higher the stimulation energy, the greater the number of muscular fibres (motor units) working. Try to gradually increase the energy level during a treatment session.
<b>Option 2+2</b>	If the patient experiences discomfort and pain during stimulation, the 2+2 function allows this programme (ch 1+2) to be combined with a pain relief programme (ch 3+4).

### BACK -TRUNK / STABILIZATION (30 MIN)

	Warm-up	Contraction	Active rest	Final recovery
Frequency	6 Hz	40 Hz	4 Hz	3 Hz
Duration of ramp-up	1.5 s	2 s	0.5 s	1.5 s
Duration of phase	2 min	6 s	12 s	3 min
Duration of ramp-down	2 s	1 s	0.5 s	3 s

<b>ATROPHY (MOD. FREQUENCY)</b>	
<b>When?</b>	Use it on weak or immobilized muscles.
<b>Why?</b>	The programme focuses on the Type I fibres to improve the stamina in the muscle.
<b>How?</b>	By modulating the frequency for a period of time that is equal to the contraction time, the contraction feels more comfortable. The contraction starts with lower freq and ends with the highest frequency.
<b>Example</b>	Increased strength for shoulder muscles by stimulation on the Rhomboid muscles. Post operative stimulation on Quadriceps, the medial vastus in particular.
<b>Electrodes</b>	Electrodes positioned as indicated according to the stimulated muscle. See suggested electrodes placements on the fold out cover.
<b>Intensity</b>	Maximum tolerable stimulation energy, which is one of the key factors for the effectiveness of the treatment. The higher the stimulation energy, the greater the number of muscular fibres (motor units) working. Try to gradually increase the energy level during a treatment session.
<b>Option 2+2</b>	If the patient experiences discomfort and pain during stimulation, the 2+2 function allows this programme (ch 1+2) to be combined with a pain relief programme (ch 3+4).

<b>ATROPHY (MOD. FREQUENCY) (30 MIN)</b>				
	Warm-up	Contraction	Active rest	Final recovery
Frequency	6 Hz	25-40 Hz	4 Hz	3 Hz
Duration of ramp-up	1.5 s	2 s	0.5 s	1.5 s
Duration of phase	2 min	4 s	8 s	3 min
Duration of ramp-down	2 s	1 s	0.5 s	3 s

**FORCE (MOD. FREQUENCY)**

<b>When?</b>	To maintain muscular strength during immobilization, or to increase muscular strength for a specific muscle/muscle group.
<b>Why?</b>	The programme focuses on the Type II fibres to improve stamina in the muscle.
<b>How?</b>	By modulating the frequency for a period of time that is equal to the contraction time, the contraction feels more comfortable.
<b>Example</b>	Increased strength for hip abductors. Stimulation on Quadriceps, in case of immobilization.
<b>Electrodes</b>	Electrodes positioned as indicated according to the stimulated muscle. See suggested electrodes placements on the fold out cover.
<b>Intensity</b>	Maximum tolerable stimulation energy, which is one of the key factors for the effectiveness of the treatment. The higher the stimulation energy, the greater the number of muscular fibres (motor units) working. Try to gradually increase the energy level during a treatment session.
<b>Option 2+2</b>	If the patient experiences discomfort and pain during stimulation, the 2+2 function allows this programme (ch 1+2) to be combined with a pain relief programme (ch 3+4).

**FORCE (MOD. FREQUENCY) (30 MIN)**

	Warm-up	Contraction	Active rest	Final recovery
Frequency	6 Hz	35-60 Hz	4 Hz	3 Hz
Duration of ramp-up	1.5 s	3 s	0.5 s	1.5 s
Duration of phase	2 min	8 s	15 s	3 min
Duration of ramp-down	2 s	1 s	0.5 s	3 s

<b>INCREASE CIRCULATION</b>	
<b>When?</b>	As a preparation before, or as recovery after, muscular activity. Could be used to increase the circulation locally in the tissue/ muscle.
<b>Why?</b>	To increase the elimination of lactic acid and waste products and to reduce stiffness in the muscles.
<b>How?</b>	This programme creates muscular twitches in order to increase the blood flow to the muscles, increase resistance strength and recover capacity.
<b>Example</b>	Put big electrodes on large muscle groups, for example Quadriceps and Hamstrings.
<b>Electrodes</b>	Electrodes positioned as indicated according to the stimulated muscle. See suggested electrodes placements on the fold out cover.
<b>Intensity</b>	The stimulation should feel pleasant and give visible muscle twitches. Visible muscle twitches are the key factor for the effectiveness of the treatment.
<b>Option 2+2</b>	If the patient experiences discomfort and pain during stimulation, the 2+2 function allows this programme (ch 1+2) to be combined with a pain relief programme (ch 3+4).

<b>INCREASE CIRCULATION (20 MIN)</b>	
<b>Frequency</b>	8 Hz continuous stimulation

## MUSCLE LESION

<b>When?</b>	If you have a more or less significant muscle lesion, this programme could be used as a complement to the usual treatment (cold, compression, etc). Mobilization of the limb must be allowed.
<b>Why?</b>	To facilitate the muscle fibres for a faster return to normal activity.
<b>How?</b>	The programme is designed to ramp-up and stimulate the muscle gradually in order to prevent stress on the muscle fibres.
<b>Example</b>	Stimulation on the back of the thigh (hamstrings).
<b>Electrodes</b>	Electrodes positioned as indicated according to the stimulated muscle. See suggested electrodes placements on the fold out cover.
<b>Intensity</b>	Maximum tolerable stimulation energy, which is one of the key factors for the effectiveness of the treatment. The higher the stimulation energy, the greater the number of muscular fibres (motor units) working. Try to gradually increase the energy level during a treatment session.
<b>Option 2+2</b>	If the patient experiences discomfort and pain during stimulation, the 2+2 function allows this programme (ch 1+2) to be combined with a pain relief programme (ch 3+4).

MUSCLE LESION (30 MIN)				
	Warm-up	Contraction	Active rest	Final recovery
Frequency	6 Hz	40 Hz	4 Hz	3 Hz
Duration of ramp-up	1.5 s	6 s	1.5 s	1.5 s
Duration of phase	2 min	3 s	10 s	3 min
Duration of ramp-down	2 s	1.5 s	1.5 s	3 s

<b>MOTOR POINT</b>	
<b>When?</b>	This programme (3 Hz continuous), combined with the use of a motor point pen, allows the user to determine the exact location of the motor point on each individual.
<b>Why?</b>	An electrode placed on the motor point ensures the most effective and comfortable stimulation.
<b>How?</b>	Connect the motor point pen to the red connector of the cable (see connections in the User manual). Connect the black connector to an electrode and place it on the skin. Put some gel on the skin area where you will search for the motor point and increase the intensity. The right spot is where you receive the strongest and most visible muscular twitch. Mark out the spot with a skin marker pencil.
<b>Intensity</b>	Very progressively increase the energy until a value between 5 and 15 is reached, while continuously moving the pen tip over the gel layer, but without ever losing contact with the skin gel, to avoid triggering off an electrode fault message.

<b>MOTOR POINT (15 MIN)</b>	
Frequency	3 Hz continuous stimulation

## 2) Pain relief

**TENS 100 HZ OR 80 HZ**

<b>When?</b>	TENS is the first choice for acute and long-term pain, both neurogenic and nociceptive. Good results are observed, above all when the pain originates in joints, the skeleton, muscles, skin, viscera or nervous system.
<b>Why?</b>	For an acute pain you need pain relief to take part in daily activities and avoid enter a vicious circle. For chronic pain you need a way to handle your daily life. TENS has no known side effects and is easy to use.
<b>How?</b>	TENS utilises the nervous system's own pain relief mechanisms and is based on the Gate control theory. This theory states that stimulating A-Beta-fibres, the sensory nerve fibres that conduct pressure, touch and vibration, inhibits impulse transfer in the pain pathways. Frequency: 80 or 100 Hz. Pulsewidth: appropriate to the sensitivity of the treated area and the patient. You can choose between 100 Hz: level 1 - 30 µs (very sensitive), level 2 - 50 µs, level 3 - 70 µs or 80 Hz: 180 µs (normal).
<b>Example</b>	Neck pain, shoulder pain, elbow pain, rheumatic pain, fractured rib, lumbago, menstrual pain, phantom limb pain, hip pain and osteoarthritic pain.
<b>Electrodes</b>	Neck pain, shoulder pain, elbow pain, rheumatic pain, fractured rib, lumbago, menstrual pain, phantom limb pain, hip pain and osteoarthritic pain.
<b>Intensity</b>	The stimulation should give a strong, but pleasant paraesthesia - tingling. Adjust the intensity level during the treatment, because this is one of the key factors for the effectiveness of the treatment.
<b>Option 2+2</b>	As a rule, the electrodes should be placed on or near the painful area, or over an area segmentally related to the painful area. See suggested electrodes placements on the fold out cover. It is important to ensure that the patient has normal sensitivity in the area where the electrodes are placed.

TENS			
Frequency	Level	Pulse width	Treatment time
100 Hz	1	30 µs	20 min
100 Hz	2	50 µs	20 min
100 Hz	3	70 µs	20 min
80 Hz	-	180 µs	30 min

## FREQUENCY MODULATED TENS

<b>When?</b>	TENS is the first choice for acute and long-term pain, both neurogenic and nociceptive. Good results are observed, above all when the pain originates in joints, the skeleton, muscles, skin, viscera or nervous system.
<b>Why?</b>	For an acute pain you need pain relief to take part in daily activities and avoid enter a vicious circle. For chronic pain you need a way to handle your daily life. TENS has no known side effects and is easy to use.
<b>How?</b>	TENS utilises the nervous system's own pain relief mechanisms and is based on the Gate control theory. This theory states that stimulating A-Beta-fibres, the sensory nerve fibres that conduct pressure, touch and vibration, inhibits impulse transfer in the pain pathways. To avoid adaption, Frequency modulated TENS have a continuous variation of the stimulation frequency. Pulsewidth: appropriate to the sensitivity of the treated area and the patient.
<b>Example</b>	Neck pain, shoulder pain, elbow pain, rheumatic pain, lumbago, menstrual pain, phantom limb pain, hip pain and osteoarthritic pain.
<b>Electrodes</b>	As a rule, the electrodes should be placed on or near the painful area, or over an area segmentally related to the painful area. See suggested electrodes placements on the fold out cover. It is important to ensure that the patient has normal sensitivity in the area where the electrodes are placed.
<b>Intensity</b>	The stimulation should give a strong, but pleasant paraesthesia - tingling. Adjust the intensity level during the treatment, because this is one of the key factors for the effectiveness of the treatment.
<b>Option 2+2</b>	If the patient experiences discomfort and pain during stimulation, the 2+2 function allows this programme (ch 1+2) to be combined with a pain relief programme (ch 3+4).

### FREQUENCY MODULATED TENS (20 MIN)

Frequency	50-150 Hz
Pulse width	50 $\mu$ s
Modulation time	2 s

## ENDORPHINIC OR BURST TENS

<b>When?</b>	The treatment Endorphinic or Burst TENS is usually most effective for radiating (projected) pain in the arms and legs (rhizopathy), for conditions with reduced or changed sensory of touch, for deep muscular pain or when the post-treatment of TENS is too short.
<b>Why?</b>	For an acute pain you need pain relief to take part in daily activities and not enter a vicious circle. For chronic pain you need a way to handle your daily life. Endorphinic or Burst TENS has no known side effects and is easy to use.
<b>How?</b>	Endorphinic or Burst TENS utilises the nervous system's own pain relief mechanisms and alleviates pain by stimulating muscles to release the body's own morphine-like substances, endorphins from hypothalamus. The muscle twitches also increase the local blood circulation. Frequency: 2-5 Hz.
<b>Example</b>	Mononeuropathy, central pain, cervical rhizopathy, sciatica and knee pain.
<b>Electrodes</b>	Place the electrodes on a muscle in the painful area so a visible twitch occurs, or on acupuncture points in the painful area. See suggested electrodes placements on the fold out cover.
<b>Intensity</b>	The stimulation should feel pleasant and give visible muscle twitches. Visible muscle twitches are the key factor for the effectiveness of the treatment.
<b>Option 2+2</b>	If the patient experiences discomfort and pain during stimulation, the 2+2 function allows this programme (ch 1+2) to be combined with a pain relief programme (ch 3+4).

### ENDORPHINIC (20 MIN)

Frequency	5 Hz continuous stimulation
Pulse width	200 $\mu$ s

### BURST TENS (20 MIN)

Frequency	2 Hz (2 pulse trains per sec. with an internal frequency of 80 Hz)
Pulse width	180 $\mu$ s

## DECONTRACTION

<b>When?</b>	To obtain a reduction in the muscular tension.
<b>Why?</b>	Studies have revealed that muscular twitches produced by a very low frequency (1 Hz) effectively reduces the muscle tension of the stimulated muscles.
<b>How?</b>	Torticollis, lumbago.
<b>Example</b>	Electrodes positioned according to the muscles where you want to reduce tension. See suggested electrodes placements on the fold out cover.
<b>Electrodes</b>	The stimulation should feel pleasant and give visible muscle twitches. Visible muscle twitches are the key factor for effectiveness of the treatment.
<b>Intensity</b>	If the patient experiences discomfort and pain during stimulation, the 2+2 function allows this programme (ch 1+2) to be combined with a pain relief programme (ch 3+4).
<b>Option 2+2</b>	Yes. If the patient experience discomfort during stimulation, the 2+2 function allows this program (ch 1+2) to be combined with a TENS program (ch 3+4) to produce the muscular twitches more quickly and more comfortably.

### DECONTRACTION (20 MIN)

Frequency	1 Hz continuous stimulation
-----------	-----------------------------

## PULSE WIDTH MODULATED TENS

<b>When?</b>	TENS is the first choice for acute and long-term pain, both neurogenic and nociceptive. Good results are observed above all when the pain originates in joints, the skeleton, muscles, skin, viscera or nervous system. Pulse width modulated TENS also gives a massage effect on muscles like the Trapezius.
<b>Why?</b>	For an acute pain you need pain relief to take part in daily activities and avoid enter a vicious circle. For chronic pain you need a way to handle your daily life. TENS has no known side effects and is easy to use.
<b>How?</b>	TENS utilises the nervous system's own pain relief mechanisms and is based on the Gate control theory. This theory states that stimulating A-Beta-fibres, the sensory nerve fibres that conduct pressure, touch and vibration, inhibits impulse transfer in the pain pathways. Pulse width modulated TENS is a type of stimulation where the pulse width varies continuously. This creates an undulating sensation which can be more pleasant than a constant pulse width.
<b>Example</b>	Neck pain, shoulder pain, trapezius pain, elbow pain, rheumatic pain, lumbago, menstrual pain, phantom limb pain, hip pain and osteoarthritic pain.
<b>Electrodes</b>	As a rule, the electrodes should be placed on or near the painful area, or over an area segmentally related to the painful area. See suggested electrodes placements on the fold out cover. It is important to ensure that the patient has normal sensitivity in the area where the electrodes are placed.
<b>Intensity</b>	The stimulation should give a strong, but pleasant paraesthesia - tingling. Adjust the intensity level during the treatment, because this is one of the key factors for the effectiveness of the treatment.
<b>Option 2+2</b>	If the patient experiences discomfort and pain during stimulation, the 2+2 function allows this programme (ch 1+2) to be combined with a pain relief programme (ch 3+4).

PULSE WIDTH MODULATED TENS (30 MIN)	
Frequency	80 Hz
Pulse width	70-180 $\mu$ s
Modulation time	2 s

## MIXED TENS

<b>When?</b>	TENS is the first choice for acute and long-term pain, both neurogenic and nociceptive. Good results are observed above all when the pain originates in joints, the skeleton, muscles, skin, viscera or nervous system. Mixed TENS is also called Han-stimulation. Stimulation frequencies switch every three seconds, giving a combination of 80 Hz and 2 Hz stimulation, which can offer a more effective treatment.
<b>Why?</b>	For an acute pain you need pain relief to take part in daily activities and avoid enter a vicious circle. For chronic pain you need a way to handle your daily life. TENS has no known side effects and is easy to use.
<b>How?</b>	Mixed TENS utilises the nervous system's own pain relief mechanisms and is based on the Gate control theory (see TENS) and also on the release of the body's own morphine-like substances, endorphins (see Burst TENS).
<b>Example</b>	Neck pain, shoulder pain, lumbago, hip pain and thigh pain.
<b>Electrodes</b>	Place the electrodes on a muscle in the painful area so that a visible contraction occurs. See suggested electrode placements on the fold out cover.
<b>Intensity</b>	The stimulation should give both a strong, but pleasant paraesthesia - tingling and as well visible muscle twitches. Note: This programme has separate intensity levels. Start to adjust the energy level for 80 Hz to a tingling sensation, then again for 2 Hz in order to have visible muscle twitches.
<b>Option 2+2</b>	If the patient experiences discomfort and pain during stimulation, the 2+2 function allows this programme (ch 1+2) to be combined with a pain relief programme (ch 3+4).

### MIXED TENS (30 MIN)

Frequency	80 Hz 3 s / 2 Hz 3 s
Pulse width	180 $\mu$ s

# IV Specific Treatment

## 1. Direct current

### A. Precautions prior to Iontophoresis treatment

**ATTENTION:** Do not perform the treatment if the patient is suffering or has suffered from asthma, hay fever, food allergy, eczema, allergy to penicillin or aspirin. Do not carry out the treatment on allergic patients, whatever form their allergy may take: hay fever, eczema, or food allergy. The more likely the medicinal product is to cause strong reactions in an allergic subject (e.g. aspirin), the more vigilant one should be.

**ATTENTION:** Make sure that the medicine is not contraindicated. Iontophoresis treatment must not be performed if the patient has a disease or is taking other treatments that are listed among the contraindications for the ionized medicine.

**ATTENTION:** Stop the treatment immediately and do not repeat it with the same medicine if a local allergic reaction is identified. Do not repeat Iontophoresis treatment if any local allergic reaction, however mild, was observed during the last treatment.

**ATTENTION:** No Iontophoresis treatment near a metal implant. Electrodes for Iontophoresis treatment must not be placed close to metallic bone or joint implants (prosthesis or bone fixing).

### B. Preparing the patient and the area to be treated by Iontophoresis

1. Thoroughly clean the area of skin to be treated, then rinse and dry.
2. Correct cleaning of the skin is not enough. It must also be degreased with a fat solvent (such as ether) applied to swabs.

**ATTENTION:** Do not shave the area of skin onto which the electrodes are placed. Hair does not interfere with Iontophoresis treatment. If treatment is done in an area where hair is shaved, there is a risk of causing small skin wounds. These wounds form points of low electrical resistance where the current will flow preferentially.

3. Place the patient in a relaxed position so that he moves as little as possible during treatment.

### C. Preparing the electrodes and solution of ionized medicine

1. Apply the solution of ionized medicine to a dry electrode previously rinsed with distilled water.
2. Apply the ionized medicinal solution to the electrode of the same polarity. In this way, the medicinal ions are repelled from that electrode and attracted to the other with the opposite polarity.
3. In order to make the circuit conductive, the active electrode has been impregnated with the solution of ionized medicine and the inactive electrode has to be soaked with a conductive substance of the therapist's choice: a conducting gel, physiological liquid, or simply tap water.

### D. Attaching the electrodes

1. Place the active electrode on the area to be treated. If the area to be treated is painful, find the chosen pain point by palpation and centre the active electrode on that point.

**ATTENTION:** Avoid placing the active electrode over scarred areas. Unless the Iontophoresis treatment is intended to soften a scar or improve a keloid, avoid placing the active electrode on an area of skin with scarring.

**ATTENTION:** Do not place electrodes over skin wounds, however slight. Apart from special forms of Iontophoresis treatment, such as antibiotic therapy for instance, only place the electrodes over healthy, intact skin with no lesions, however slight.

2. When attaching the electrodes it is important to make sure that their entire surface area is applied to the skin. Just applying a strap passing through the centre of the electrode and leaving the outer edges unattached is inadmissible. Use the widest possible strap, use several straps or several turns of the same strap or even use adhesive tape to fix the sides of the electrodes properly.

**ATTENTION:** Make sure there is never any contact between a metal component and the skin. If the connector of an electrode comes into contact with the skin, the current will flow preferentially through that point of low impedance. As this contact has a very small surface area, the density of the electricity will be very high, resulting in an electric burn.

3. If possible, place the inactive electrode at right angles to the active electrode. There has been no study on how the positioning of the two electrodes in relation to each other influences the efficacy of Iontophoresis treatment. However, the depth of penetration should logically be greater if the direction of the electric field is perpendicular to the surface of the skin rather than oblique or longitudinal.

#### E. During treatment

**ATTENTION:** Do not move or remove the electrodes without stopping the treatment first. Physio 5 is programmed so that the current increases gradually at the start of the treatment and decreases gradually at the end or when the treatment is stopped. This means that there can be no excitation phenomenon and the patient will never be surprised by a shock or a painful electrical discharge. If, by contrast, the electrodes are disconnected the sudden break in the circuit may give rise to an excitation phenomenon.

1. Ask the patient to move as little as possible during the treatment and not to remove the electrodes. For the same reasons as in the previous point.
2. Warn the patient that a pricking sensation from the electrodes is normal and harmless. This is a normal effect of the galvanic current which has nothing to do with burning.
3. If there is an electrode fault during treatment. The Physio 5 measures the impedance of the circuit and, when this is too high, the equipment stops and indicates "ELECTRODE FAULT" as well as the number of the channel on which there is a problem. There are a number of possible reasons for this safety and efficacy check system coming into operation:
  - electrode disconnected
  - channel reversal
  - defective electrode
  - poor connection
  - defective cable
  - solution not conducting (non-ionizable medicine or concentration too low)

#### F. After treatment

1. Thoroughly clean the skin over the treated area using tap water. During Iontophoresis treatment, acids and bases form on the electrodes and hence come into contact with the skin. If the concentration of these substances is too high and they stay on the skin for too long, chemical burns may result. It is advisable to clean the patient's skin immediately after the treatment to remove these chemical substances.
2. Clean the electrodes thoroughly with tap water, then rinse with distilled water before leaving them to dry.

## 2. Denervated

### A. Introduction

In the current state of knowledge there is nothing to indicate that electrostimulation is capable of influencing the re-innervation process of a partly or fully denervated muscle. Electrostimulation of denervated muscle fibres, however, is essential insofar as it is the only really effective means of retaining a certain trophicity and limiting the sclerosis phenomenon of these fibres throughout the duration of their possible re-innervation period. Indeed, after many months of being patient, nothing is more frustrating than to find functional trouble caused by muscles that are certainly re-innervated but with a sclerosis condition that prevents them from being used satisfactorily. If stimulation enables the amyotrophy to be limited and sclerosis of the denervated muscle to be avoided during its re-innervation period, it then becomes pointless if there is any hope of re-innervation for the denervated fibres.

The choice of form and parameters of the electrical current depend on state of denervation of the muscle: is it completely or partly denervated? Therefore, before undertaking any electrostimulation treatment on a denervated muscle, the following two questions should be answered:

1. Is there any hope of re-innervation? In other words, have the re-innervation times elapsed or not?
2. Is the muscle completely or partly denervated?

### B. Factors guiding the therapeutic approach

#### 1. Are we within the re-innervation times?

To be able to answer this question, it is essential to have the following three pieces of information:

- A The date of the injury,
- B The degree of the injury,
- C The rate of nerve fibre regeneration.
  - Interviewing the patient usually establishes how old the injury is and where it is located.
  - The rate of regeneration of an injured nerve is approximately 1 millimetre per day, i.e. 3 centimetres a month.
  - The following elementary calculation gives the re-innervation times:

$$\frac{\text{Distance in cm between nerve fibre injury and the motor point of the denervated muscle}}{\text{Rate of nerve fibre regeneration (} = 3 \text{ cm per month)}} = \text{Re-innervation time}$$

#### 2. Total or partial denervation?

How can we find out if the muscle is partly or totally denervated?

- An electromyogram examination is of course preferable but it must be recent and the results must be passed on to the physiotherapist, which does not always happen in day to day practice.
- Muscular testing is often worthwhile. However, with certain muscles, especially if there are only very few innervated fibres left, the really analytical contraction of the muscle is difficult to obtain because of the inevitable activity of the agonist muscles.
- Nevertheless, there is a simple and easily reproducible way to find out the state of denervation of a muscle. Biphase rectangular micro-pulses (lasting between 0.15 and 0.35 ms) are only capable of exciting the nerves but not of directly exciting the denervated muscle fibres. It is sufficient, therefore, to test by


means of a disuse atrophy treatment. If no response is observed in spite of significant current strengths, the muscle can then be considered as completely denervated; if, on the other hand, a contraction, even of low intensity, is achieved, then the muscle is partly denervated.

### C. Practical therapeutic approach

It is therefore actually easy to find out the two fundamental factors that will guide our therapeutic approach:

- There is hope of re-innervation or, on the contrary, denervation is final.
- The muscle is partly or totally denervated.

Four situations can thus arise:

The practical therapeutic approach must be adapted to each situation:

#### Situation 1: Total denervation outside the time

Electrostimulation by means of the Denervated programmes is pointless, since a muscle definitively without any innervation will always end up atrophying and sclerosing.

#### Situation 2: Partial denervation outside the time

It is not possible to avoid atrophy and sclerosis of muscle fibres that are definitively denervated. Stimulation of these fibres by means of the Denervated programmes is therefore not indicated here. It is possible, however, to work on the innervated part of the muscle, by means of neurostimulation rectangular biphasic micropulses in order to achieve compensatory hypertrophy of the innervated fibres.

#### Situation 3: Total denervation within the time

Pending possible re-innervation, it is important to prevent atrophy as much as possible and limit the sclerosis phenomenon. Stimulation of muscles deprived of innervation, by means of wide rectangular pulses in the Denervated programmes is the preferred technique here.

#### Situation 4: Partial denervation within the time

It is important to try and prevent atrophy and to limit the phenomenon of sclerosis of the denervated fibres; to do this it is necessary to use the triangular gradient pulses in the Denervated programmes. Depending on the circumstances it may also be worthwhile working on the innervated part of the muscle using the rectangular biphasic micro-pulses in the neurostimulation programmes.

### 3) Aesthetic

TONING	
<b>Indication</b>	To be used initially in order to tone up and prepare the muscles before more intensive firming work. This type of training session is of moderate intensity, making it suitable for everyone, and does not cause muscle soreness.
<b>Effects</b>	To tone the muscles.
<b>Electrodes</b>	Electrodes positioned as indicated according to the stimulated muscle. See suggested electrodes placements on the fold out cover.
<b>Intensity</b>	Maximum tolerable stimulation energy, which is one of the key factors for the effectiveness of the treatment. The higher the stimulation energy, the greater the number of muscular fibres (motor units) working. Try to gradually increase the energy level during a treatment session.
<b>Option 2+2</b>	If the patient experiences discomfort and pain during stimulation, the 2+2 function allows this programme (ch 1+2) to be combined with a pain relief programme (ch 3+4).

TONING (25 MIN)				
	Warm-up	Contraction	Active rest	Final recovery
Frequency	6 Hz	18 Hz	4 Hz	3 Hz
Duration of ramp-up	1.5 s	1.5 s	0.5 s	1.5 s
Duration of phase	2 min	6 s	10 s	3 min
Duration of ramp-down	2 s	0.75 s	0.5 s	3 s

<b>FIRMING</b>	
<b>Indication</b>	To be used as the main treatment for muscle firming. Firming programmes represent exercise of medium intensity, and are aimed to improve muscular density without causing muscle soreness.
<b>Effects</b>	To regain muscle firmness and restore the support function of the muscles.
<b>Electrodes</b>	Electrodes positioned as indicated according to the stimulated muscle. See suggested electrodes placements on the fold out cover.
<b>Intensity</b>	Maximum tolerable stimulation energy, which is one of the key factors for the effectiveness of the treatment. The higher the stimulation energy, the greater the number of muscular fibres (motor units) working. Try to gradually increase the energy level during a treatment session.
<b>Option 2+2</b>	If the patient experiences discomfort and pain during stimulation, the 2+2 function allows this programme (ch 1+2) to be combined with a pain relief programme (ch 3+4).

<b>FIRMING (22 MIN)</b>				
	Warm-up	Contraction	Active rest	Final recovery
Frequency	6 Hz	22 Hz	4 Hz	3 Hz
Duration of ramp-up	1.5 s	1.5 s	1 s	1.5 s
Duration of phase	2 min	8 s	5 s	3 min
Duration of ramp-down	2 s	0.75 s	1 s	3 s

## 4) Incontinence

**STRESS INCONTINENCE**

<b>Indication</b>	The sphincter urethra is deficient and cannot remain closed in the event of a sudden and significant increase in abdominal pressure (exertion, coughing, etc.).
<b>Effects</b>	The aim of this programme is to strengthen the sphincter muscle of the bladder. The objective is therefore to produce tetanic contractions of the paraurethral components of the striated muscle of the pelvic floor, using the optimum tetanization frequencies of the fast fibres.
<b>Electrodes</b>	Use of an intravaginal probe.
<b>Intensity</b>	Use of the maximum energy tolerated by the patient in order to obtain the greatest possible spatial recruitment. The current intensity is increased regularly every 3 or 4 contractions throughout the session. The therapist plays a decisive role in reassuring the patient and obliging her to work with the most powerful possible contractions.
<b>Option 2+2</b>	No

**STRESS INCONTINENCE (20 MIN)**

	Contraction	Active rest
Frequency	75 Hz	0 Hz
Duration of ramp-up	1.5 s	0 s
Duration of phase	4 s	12 s
Duration of ramp-down	1.5 s	0 s

## URGE INCONTINENCE

<b>Indication</b>	The bladder contracts abnormally (detrusor overactivity) and presses on the urine, increasing the pressure within the bladder.
<b>Effects</b>	This treatment relies on the reduction of detrusor activity by stimulation of an inhibitory reflex from sensory nerve endings in the perineal region. The electrical parameters must therefore be created so as to excite these myelinated afferent nerve fibres at the frequency that produces optimum activation of the inhibitory reflex.
<b>Electrodes</b>	Use of an intravaginal probe.
<b>Intensity</b>	Gradually increase the energies until the patient feels the stimulation. Then increase the energies again to a value equal to three times that of the perception threshold.
<b>Option 2+2</b>	No

### URGE INCONTINENCE (30 MIN)

<b>Frequency</b>	5 Hz continuous stimulation
<b>Pulse width</b>	150 $\mu$ s

## MIXED INCONTINENCE

<b>Indication</b>	Combination of urge and stress incontinence in greater or lesser proportions.
<b>Effects</b>	This programme treats both aspects of this form of incontinence at the same time. Firstly, using tetanic contractions at the frequency of fast fibres (75 Hz), it strengthens the paraurethral components of the striated muscle of the pelvic floor, so increasing the pressure of urethral closure. Secondly, during the resting phases between contractions, it inhibits the activity of the smooth muscle of the bladder using very low frequencies (5 Hz).
<b>Electrodes</b>	Use of an intravaginal probe.
<b>Intensity</b>	Use of the maximum energy level tolerated during the tetanic contraction phases to obtain the maximum possible space recruitment and therefore maximum possible efficacy. The intensity will be increased regularly during the session, every 3 or 4 contractions. During the rest phase, the low-frequency intensity should be adjusted to at least three times the intensity of the perception threshold.
<b>Option 2+2</b>	No

MIXED INCONTINENCE (30 MIN)		
	Contraction	Active rest
Frequency	75 Hz	5 Hz
Duration of ramp-up	1.5 s	0.5 s
Duration of phase	4 s	23 s
Duration of ramp-down	0.75 s	0.5 s

## POST PARTUM PREVENTION

<b>Indication</b>	Labour causes considerable trauma to the pelvic region. The consequences of this trauma are various: strained muscle, torn muscle, partial denervation, loss of body image, loss of strength and control of the striated muscles of the pelvic floor, etc.
<b>Effects</b>	Incontinence is a relatively common result of this situation, which is why prophylactic pelvic re-training treatment by neuromuscular electrostimulation is indicated.
<b>Electrodes</b>	Use of an intravaginal probe.
<b>Intensity</b>	Use of the maximum electrical intensity tolerated by the patient in order to achieve the greatest possible space recruitment. The electrical intensity is increased regularly every 3 or 4 contractions throughout the session.
<b>Option 2+2</b>	No

### POST PARTUM PREVENTION (20 MIN)

	Contraction	Active rest
Frequency	50 Hz	0 Hz
Duration of ramp-up	1.5 s	0 s
Duration of phase	5 s	10 s
Duration of ramp-down	0.75 s	0 s

# V Conditioning

## STRENGTH

<b>Indication</b>	For competitive athletes who practice a discipline that requires strength and speed. Recommended use is three times a week in order to improve the power.
<b>Effects</b>	This programme increases the maximum strength and the rate of muscular contractions. Integrate this programme into your voluntary training programme.
<b>Electrodes</b>	Electrodes positioned as indicated according to the stimulated muscle. See suggested electrodes placements on the fold out cover.
<b>Intensity</b>	Maximum tolerable stimulation energy, which is one of the key factors for the effectiveness of the treatment. The higher the stimulation energy, the greater the number of muscular fibres (motor units) working. Try to gradually increase the energy level during a treatment session.
<b>Option 2+2</b>	If the patient experiences discomfort and pain during stimulation, the 2+2 function allows this programme (ch 1+2) to be combined with a pain relief programme (ch 3+4).

### STRENGTH (33 MIN)

	Warm-up	Contraction	Active rest	Final recovery
Frequency	5 Hz	75 Hz	4 Hz	3 Hz
Duration of ramp-up	1.5 s	1.5 s	0.5 s	1.5 s
Duration of phase	5 min	4 s	19 s	10 min
Duration of ramp-down	2 s	0.75 s	0.5 s	3 s

## EXPLOSIVE STRENGTH

<b>Indication</b>	For athletes who practice a discipline in which explosive strength is an important factor in performance. To increase the capacity for instantaneous maximum force. Recommended use is three times a week.
<b>Effects</b>	This programme increases the speed at which the level of strength is attained and improves the efficiency of explosive actions like jumping, sprinting etc.
<b>Electrodes</b>	Electrodes positioned as indicated according to the stimulated muscle. See suggested electrodes placements on the fold out cover.
<b>Intensity</b>	Maximum tolerable stimulation energy, which is one of the key factors for the effectiveness of the treatment. The higher the stimulation energy, the greater the number of muscular fibres (motor units) working. Try to gradually increase the energy level during a treatment session.
<b>Option 2+2</b>	If the patient experiences discomfort and pain during stimulation, the 2+2 function allows this programme (ch 1+2) to be combined with a pain relief programme (ch 3+4).

### EXPLOSIVE STRENGTH (32 MIN)

	Warm-up	Contraction	Active rest	Final recovery
Frequency	5 Hz	104 Hz	1 Hz	3 Hz
Duration of ramp-up	1.5 s	0.75 s	0.5 s	1.5 s
Duration of phase	5 min	3 s	28 s	10 min
Duration of ramp-down	2 s	0.5 s	0.5 s	3 s

## HYPERTROPHY

<b>Indication</b>	For body building enthusiasts and athletes who want to increase their muscle mass. Combine this programme with the voluntary training programme. Recommended use is three to five sessions a week.
<b>Effects</b>	This programme increases the volume of the stimulated muscles and improves the muscle resistance.
<b>Electrodes</b>	Electrodes positioned as indicated according to the stimulated muscle. See suggested electrodes placements on the fold out cover.
<b>Intensity</b>	Maximum tolerable stimulation energy, which is one of the key factors for the effectiveness of the treatment. The higher the stimulation energy, the greater the number of muscular fibres (motor units) working. Try to gradually increase the energy level during a treatment session.
<b>Option 2+2</b>	If the patient experiences discomfort and pain during stimulation, the 2+2 function allows this programme (ch 1+2) to be combined with a pain relief programme (ch 3+4).

### HYPERTROPHY (31 MIN)

	Warm-up	Contraction	Active rest	Final recovery
Frequency	5 Hz	45 Hz	8 Hz	3 Hz
Duration of ramp-up	1.5 s	1.5 s	5 s	1.5 s
Duration of phase	5 min	4 s	8 s	10 min
Duration of ramp-down	2 s	1 s	0 s	3 s

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© 2010 DJO, LLC
Cefar-Compex® is a registered trademark of DJO, LLC.
DJO, LLC
Vista, CA 92081 U.S.A.
MDSS GmbH