



# ARTROMOT®-S3 COMFORT

for the shoulder joint

## Technical Data:

### Models:

ARTROMOT®-S3

ARTROMOT®-S3 COMFORT

### Range of motion:

Adduction/abduction	0° / 30° / 175°
Internal/external rotation	90° / 0° / 90°
Elevation	0° / 30° / 175°
Horizontal adduction/abduction (manual)	0° / 0° / 120°

### Functions of ARTROMOT®-S3:

Motor control	ON/OFF for isolated motion
Synchronized motors	ON/OFF
Speed (ABD/ADD)	1% – 100%
Timer	1 – 300 minutes / continuous operation
Pauses	0 – 30 seconds
Reverse on load in steps	1 – 25 / patient safety feature
Total therapy time	added sum of the therapy sessions
Warm up protocol	for gradual progression to the maximum range of motion
Patient chip card	for storage of the therapy parameters
New patient	activates the default settings for new patients
Service menu	

### Additional functions of ARTROMOT®-S3 COMFORT:

Headrest	for comfortable, optimal patient posture
Shoulder positioning	prevents evading movements
Oscillation	higher rate of motion in the EROM
Stretching (ABD or I-ROT or E-ROT)	careful and smooth increase of the range of motion
ISO ABD ROT protocol	for separate exercising in both directions

### Physical:

Weight:	25 kg (55.1 lb)
Dimensions	87.5 x 57.5 x 58 cm (34.4 x 22.6 x 22.8 in)
Transport dimensions	37 x 57,5 x 58 cm (14.5 x 22.6 x 22.8 in)

### Article numbers:

ARTROMOT®-S3	80.00.023
ARTROMOT®-S3 COMFORT	80.00.024
ARTROMOT® chip card	0.0034.048